MALDI-TOF/TOF Self User Training

2013. 12. 31

Jinyoung Kim

UNIST Central Research Facilities (UCRF)

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3. Data Upload & Download System (Core FTP)

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Data Upload & Download System

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Site Manager dialog box for external PC connection. The 'Host / IP / URL' field is highlighted with a red box and contains '10.24.9.42'. The 'Port' field is also highlighted with a red box and contains '22'. Other fields include Site Name (New_Site), Username (sophia), Password (masked), Timeout (60), Retries (66), and Connection (FTP).

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Site Manager dialog box for UCRF internal PC connection. The 'Host / IP / URL' field is highlighted with a red box and contains '100.100.100.100'. The 'Port' field is also highlighted with a red box and contains '22'. Other fields include Site Name (New_Site), Username (sophia), Password (masked), Timeout (60), Retries (66), and Connection (FTP).

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
Core FTP LE - 10.24.9.42:22

File View Sites Manage Help

250 CWD command successful. "/" is current folder.
 PWD
 257 "/" is current directory.
 PASV
 227 Entering Passive Mode (10,24,9,42,146,187).
 LIST
 Connect socket #1020 to 10.24.9.42, port 37563...
 150 Opening data conn
 226 File sent ok
 Transferred 4,523 byte: 분석 PC or My PC

Server PC

Filename	Size	Date
..		09/13/13 09:07
Fax		08/13/13 11:28
Malvern Instruments		09/19/12 10:43
SAP		10/03/13 17:24
Scanned Documents		08/13/13 11:28
경남은행		08/07/13 15:49
Self 교육.xlsx	14 KB	09/26/12 17:59
XRF_data_보정.xls	112 KB	12/27/12 13:22
거래내역조회_2009_출력.pdf	45 KB	08/07/13 16:13
동위원소실가동률계산_박지혜.xlsx	32 KB	04/26/13 18:00

Filename	Size	Date
<..>		
..		
Company		10/04/13 11:...
Labs		10/06/13 20:...
Program		09/02/13 12:...
toengineer		01/30/13 10:...
touser		10/04/13 15:...
UCRF		10/02/13 17:...
UCRF 로그복관리		10/01/13 14:...
AFM 스캐너수리 요구내역.xlsx	0	04/27/13 15:...

A Record of Equipment Use

Approval	Manager

Equipment			
University or Company			
Client	Department <small>(Supervisor)</small>		Position
	Name		Contact number
	E-mail	<small>(Please write the contact number and e-mail address in case of an outside client)</small>	
User	Department <small>(Supervisor)</small>		Position
	Name		Contact number
Date	From () Month - () Day - () Year to () Month - () Day - () Year (Total use time : hours , [: - :] minutes		
Sample		No. of samples	
Details about use of equipment			
Details about use of hazardous materials			
Equipment state (problem and repair)	Cleanliness	Good	Poor

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
What is MALDI

Matrix-assisted laser desorption/ionization (MALDI) is a soft ionization technique used in mass spectrometry, allowing the analysis of biomolecules (biopolymers such as DNA, proteins, peptides and sugars) and large organic molecules (such as polymers, dendrimers and other macromolecules), which tend to be fragile and fragment when ionized by more conventional ionization methods. It is similar in character to electrospray ionization in that both techniques are relatively soft ways of obtaining large ions in the gas phase, though MALDI produces many fewer multiply charged ions.

MALDI is a two step process. First, desorption is triggered by a UV laser beam. Matrix material heavily absorbs UV laser light, leading to the ablation of upper layer (~micron) of the matrix material. The hot plume produced during ablation contains many species: neutral and ionized matrix molecules, protonated and deprotonated matrix molecules, matrix clusters and nanodroplets. **Second, the analyte molecules are ionized (more accurately protonated or deprotonated) in the hot plume.** Ablated species may participate in the ionization of analyte, though the mechanism of MALDI is still debated.

Matrix Properties*:

- 1) Mix well with analyte, and separate molecules of analyte.
- 2) No chemical reaction and no adduct formation with analyte
- 3) Absorbption of UV
- 4) Desorption upon laser irradiation
- 5) Support ionization


- Michael Karas, Doris Bachmann and Franz Hillenkamp; Anal. Chem. 1985, 57, 2935-2939.
- Ronald C. Bavis and Brian Chait; Rapid Com. In Mass Spectr. Vol.3, No. 7, 1989
- Michael Karas, Matthias Glückmann and Jürgen Schäfer; J. Mass Spectr. 35, 1-12, (2000)


What is TOF MS

TOF MS is the abbreviation for Time of Flight Mass Spectrometry.

Charged ions of various sizes are generated on the sample slide, as shown in the diagram. A potential difference V_0 between the sample slide and ground attracts the ions in the direction shown in the diagram. The velocity of the attracted ions v is determined by the law of conservation of energy. As the potential difference V_0 is constant with respect to all ions, ions with smaller m/z value (lighter ions) and more highly charged ions move faster through the drift space until they reach the detector. Consequently, the time of ion flight differs according to the mass-to-charge ratio (m/z) value of the ion. The method of mass spectrometry that exploits this phenomenon is called Time of Flight Mass Spectrometry.

- The reflector or ion mirror compensates for the initial energy spread of ions of the same mass coming from the ion source, and improves resolution.
- Increasing of resolution and signal-to noise by focusing the ions in time and space!

Features of MALDI-TOF

- Soft ionization - analyze intact biomolecules and synthetic polymers
- Broad mass range - analyze a wide variety of biomolecules
- Simple mixtures are okay
- Relatively tolerant of buffers and salts
- Fast data acquisition
- Easy to use and maintain, no water or gas hook ups required
- High sensitivity, superior mass resolution and accuracy

Important parameters

- correct mixture using suitable solvent, matrix and cationization agent

Matrices

DHB

HCCA


Ferulic acid

IAA

Dithranol

DCTB

Hydrophilic

Hydrophobic

Polymers

Polyethylene glycol

Polydimethylsiloxane

Polytetramethylene glycol


Polymethylmethacrylate

Polystyrene







Polyisoprene


Matrix Guide to Sample Preparation

* Please refer to the attached file(MatrixSelection.pdf) for details

Matrix	Abbr.	Formular	CAS No.	Solvent	Compound
2,5-Dihydroxybenzoic acid MW = 154.12	DHB		490-79-9	THF Aceton DMF Acetonitrile Water	Peptides, Proteins Oligosaccharides Lipids Synth. Polymers Small. Org. molecules
α -Cyano-4-hydroxy cinnamic acid MW = 189.17	CCA		28166-41-8	THF DMF EtOH Acetonitrile Water	Peptides Proteins Oligosaccharides
<i>trans</i> -3,5-dimethoxy-4-hydroxy-cinnamic acid MW = 224.21	SA (Sinapinic acid)		530-59-6	DMF EtOH Acetonitrile Water	Proteins
3-Hydroxypicolinic acid MW = 139.11	3HPA		874-24-8	Acetonitrile Water	Nucleotides
1,8,9-Trihydroxy-antracene MW = 226.23	Dithranol		1143-38-0	THF Chloroform	Unpolar Synth. polymers
<i>trans</i> -Indole-3-acrylic acid MW = 187.19	IAA		530-59-6	THF Aceton	Unpolar Synth. polymers

Preferred application: fast and homogenous (flat)

- no separation
- high shot-to-shot reproducibility


Linear or Reflector / Positive or Negative


1. Reflector mode is used in the mass range up to 7 kDa, where the isotopes of the oligomers can be seen
2. Sometimes it makes sense to use the reflector mode to obtain oligomeric resolution also in the higher mass range
3. Linear mode is used for masses above 7 kDa
4. Normally polymers are observed as cation-attached species
5. Anions are only observed in the case of present negative charges in the polymer such as Poly(styrenesulfonate). In this case all cations have to be exchanged against ammonium

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
Linear / Positive

Linear / Negative

The image displays two panels from a software interface. The left panel is the 'Detection' tab, showing a 'Mass Range' slider from 0 to 10000 Da, currently set to 0-4000 Da. Below it are 'Detector Gain' and 'Sample Rate and Digitizer Settings' sections. The right panel is the 'Spectrometer' tab, showing 'High Voltage' settings and a 'Matrix Suppression' section where 'Deflection' is selected and 'Suppress up to:' is set to 500 Da. Both the 'Mass Range' and 'Matrix Suppression' sections are highlighted with red boxes.

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