# Sputtering system(DC & RF)

---

2013. 12. 27

Hyung Il Kim

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



| NO | Name                        | Description                                                                                                                                                                                                                                                                             |
|----|-----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1  | Signal Tower                | Normal, Caution, Warning Signal Lamp                                                                                                                                                                                                                                                    |
| 2  | Process Chamber Module(PM)  | Table rotation, Gun shutter, Source, View port LED, etc                                                                                                                                                                                                                                 |
| 3  | Loadlock Chamber Module(LL) | Robot                                                                                                                                                                                                                                                                                   |
| 4  | Frame module                | RF & DC Generator, MPD , EM0, Interface, Computer, etc                                                                                                                                                                                                                                  |
| 5  | Facility Requirement        | <ul style="list-style-type: none"> <li>- Main Power : 3<math>\phi</math> 220V, 50A</li> <li>- Water : 50 ~ 90 Psi, 18 to 27<math>^{\circ}</math>C, 7.5 LPM ~ 15LPM</li> <li>- Air : 75 ~ 85 Psi</li> <li>- N<sub>2</sub> : 20 ~ 35 Psi</li> <li>- Ground : Less than 0.3 ohm</li> </ul> |

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**Sorona Inc.**

**Login User**  
User: system  
Group: SYSTEM

**Module Status**  
PM: READY  
LL: READY

**Soft EMO**

**Monday 30 November 2009 20:00:29**

**Gun Selection**  
-<RF1001\_1>-  
GUN 1  
GUN 2  
-<RF1001\_2>-  
GUN 3  
GUN 4

**Talbe Rotation**  
Status: STOP  
Position: Gun1  
ROTATE STOP MOVE STOP

**Gun Shutter**  
Status: STOP  
Position: Gun1  
MOVE STOP

**PSD301W-1 [DC Power Control]**  
Set Actual  
Power: 0 0 W  
Voltage: 0.0 0.0 V  
Current: 0.0 0.0 A  
On Time: 0 0  
DC On DC Off

**PSD301W-2 [DC Power Control]**  
Set Actual  
Power: 0 0 W  
Voltage: 0.0 0.0 V  
Current: 0.0 0.0 A  
On Time: 0 0  
DC On DC Off

**R301 [Bias RF Power Control - 300W]**  
Set Actual  
Power: 0 0 W  
On Time: 0 0 sec  
Reflect: 0 0 W  
RF On RF Off

**Full Auto Control**  
Recipe Name:   
Step Name: READY  
Process Time: 0.0 / 300.0 sec  
Step Time: 0.0 / 60.0 sec  
LOAD START STOP  
Process Option:  
[x] Substrate Pre Cleaning  
[x] Auto UnLoading  
[x] Auto Vent LL

**System Interloc**  
Air  
Water  
+DC24V

**Sequence Status**  
READY  
Command  
00:00 39"  
0.0 / 1.0 sec  
Stop

**Process Timer**  
Start: 00:00:00  
End: 00:00:00  
Total: 00:00:00

**Substrate Status**  
Size: 6"  
Position: LL

**Robot Control**  
READY  
Loading  
Unloading

**LL Pumping**  
Pumping  
Venting


**PM Pumping**  
Pumping  
Venting

**Service Process** **Service Pumping** **Recipe Editor** **History** **System Config** **User Management** **System I/O Manager** **Login/out Shutdown**

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
| Window name           | Function                        | Description                             |
|-----------------------|---------------------------------|-----------------------------------------|
| Recipe control        | Recipe Select                   | Full auto recipe select                 |
|                       | Auto processing control         | Full auto Start, Pause, Stop control    |
|                       | Step name                       | Recipe Step name                        |
|                       | Process Time                    | Recipe Process Full Time                |
|                       | Step time                       | Recipe step time                        |
| DC & RF Power Control | DC & RF Power Set, Power ON/OFF | DC & RF power Set, DC On/Off, RF On/Off |
| Auto Process Option   | Auto Process End Vented         | Venting After Auto Process End          |
| Auto process message  | Auto Status                     | Auto processing status message          |

# Equipment operation

## 2) Service Pumping: PM Vacuum and Moving Control


| Window name          | Function                | Description                          |
|----------------------|-------------------------|--------------------------------------|
| PM Pumping & Venting | PM Vacuum Control       | PM Pumping, Venting control          |
| LL Robot Control     | Robot Manual Control    | Robot Retract, Extend manual control |
|                      | Robot Home              | Robot Home Position                  |
|                      | Stop                    | Robot Moving Stop                    |
| LL Pumping & Venting | Loadlock Vacuum control | Loadlock Pumping, Venting control    |
| Gun Shutter Control  | Gun Shutter Moving      | Gun Shutter Position Select & Moving |

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| Window name        | Function                                                | Description                                                                                |
|--------------------|---------------------------------------------------------|--------------------------------------------------------------------------------------------|
| Module Recipe List | Recipe Add, Delete, Edit                                | Recipe Category, Recipe Add, Recipe Delete, Recipe Edit                                    |
| Recipe Step Data   | Recipe Step Data Insert, Delete, Recipe Save, Step Move | Recipe Step Data Set, Step Insert, Step Delete, Recipe Save, Step Data Left & Right Moving |

# Equipment operation

## 4) History screen : Process & System Log, Alarm History

**Sorona Inc.**

Login User: User: system, Group: SYSTEM

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| Item                         | Name                                            | Description                                                                                                                                                              |
|------------------------------|-------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| LL Module                    | LL Module Parameter Settings                    | LL Base Pressure, LL Roughing Pressure, Robot Moving Parameter Settings, Gun shutter Moving command Settings, Pumping Time Out, Venting Time Out                         |
| PM Module                    | PM Module Parameter Settings                    | PM Base Pressure, PM Roughing Pressure, MFC Name & Range Settings, Moving Parts Time Out settings, Cryo Pumping Parameter Settings, Pumping & Venting Time Out Settings, |
| Robot                        | Robot Motor Parameter Settings                  | Robot Sensor Input/Output Channel Settings, Moving Time Settings, Moving Speed Settings, etc                                                                             |
| Gun Shutter                  | Gun Shutter Parameter Settings                  | Gun Shutter Sensor Input/Output Channel Settings, Moving Time Settings, Moving Speed Settings, etc                                                                       |
| Target Material              | Target Material List                            | Target Material Add, Edit, Delete, Update Function                                                                                                                       |
| Gun shutter & Target Mapping | Gun shutter settings & Generator/Target Mapping | Gun No – Target material - Power generator Mapping                                                                                                                       |
| RF Bias                      | R301 Bias Settings                              | RF Bias Power Settings, Pressure, On Time, Ar Flow, Stabilization Setting                                                                                                |

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**Sorona Inc.**

**Login User**  
User: system Group: SYSTEM

**Module Status**  
PM: READY LL: READY

Soft EMO

Tuesday 1 December 2009  
10:32:28

**System Manager**

| Name  | ID      | Description               |
|-------|---------|---------------------------|
| UNIST | manager | system manager (PW:UNIST) |

**Power User**

| Name  | ID    | Description            |
|-------|-------|------------------------|
| POWER | POWER | Power User (PW: POWER) |

**Operator**

| Name | ID   | Description         |
|------|------|---------------------|
| OPER | OPER | Operator (PW: OPER) |

**System Interloc**  
Air Water +DC24V

**Sequence Status**  
READY  
Command  
00:00:12"  
0.0/1.0 sec  
Stop

**Process Timer**  
Start: 00:00:00  
End: 00:00:00  
Total: 00:00:00"

**Substrate Status**  
Size: 6" Position: LL

**Robot Control**  
READY  
Loading Unloading

**LL Pumping**  
Pumping Venting

**PM Pumping**  
Pumping Venting

Service Process Service Pumping Recipe Editor History System Config User Management System I/O Manager Login/Out Shutdown

| Item           | Name             | Description       |
|----------------|------------------|-------------------|
| System Manager | Manager Group    | Add, Edit, Delete |
| Power User     | Power User Group | Add, Edit, Delete |
| Operator       | Operator Group   | Add, Edit, Delete |

# Equipment operation

## 7) System I/O Management screen : Digital Input/Output Check, Driver Status Check

| J  | Define                      | TYPE  | Ch  | Status                | Actuator          |
|----|-----------------------------|-------|-----|-----------------------|-------------------|
| 1  | R SYS DC24P INT             | DI 0  | P56 | RELEASE/DC24P         | Relay             |
| 2  | SPARE                       | DI 1  | P57 |                       |                   |
| 3  | SPARE                       | DI 2  | P58 |                       |                   |
| 4  | R SYS AIR INT               | DI 3  | P59 | RELEASE/AIR           | Sensor            |
| 5  | R SYS WATER INT             | DI 4  | P60 | RELEASE/WATER         | Sensor            |
| 6  | R PM ATM                    | DI 5  | P61 | ATM/VAC               | Sensor            |
| 7  | R PM PIN UP                 | DI 6  | P62 | UP/UNKNOWN            | Sensor            |
| 8  | R PM PIN DW                 | DI 7  | P63 | DOWN/UNKNOWN          | Sensor            |
| 9  | R PM TABLE UP               | DI 8  | P64 | UP/UNKNOWN            | Sensor            |
| 10 | R PM TABLE DOWN             | DI 9  | P65 | DOWN/UNKNOWN          | Sensor            |
| 11 | R PM GUN SHUTTER UP         | DI 10 | P66 | UP/UNKNOWN            | Sensor            |
| 12 | R PM GUN SHUTTER DOWN       | DI 11 | P67 | DOWN/UNKNOWN          | Sensor            |
| 13 | R PM GUN SHUTTER MOTOR HOME | DI 12 | P68 | HOME/UNKONW           | Sensor            |
| 14 | R PM SUB SHUTTER OPEN       | DI 13 | P69 | OPEN/UNKNOWN          | Sensor            |
| 15 | R PM SUB SHUTTER CLOSE      | DI 14 | P70 | CLOSE/UNKNOWN         | Sensor            |
| 16 | R PM RF BIAS CONNECTION     | DI 15 | P71 | CONNECTION/UNKNOWN    | Sensor            |
| 17 | R PM RF BIAS DISCONNECTION  | DI 16 | P72 | DISCONNECTION/UNKNOWN | Sensor            |
| 18 | R LL ATM                    | DI 17 | P73 | ATM/VAC               | Sensor            |
| 19 | R LL DOOR CLOSE             | DI 18 | P74 | CLOSE/OPEN            | Sensor            |
| 20 | R LL GATE OPEN              | DI 19 | P75 | OPEN/UNKNOWN          | Sensor            |
| 21 | R LL GATE CLOSE             | DI 20 | P76 | CLOSE/UNKNOWN         | Sensor            |
| 22 | R LL ROBOT LIMIT            | DI 21 | P77 | LIMEIT/UNKNOWN        | Sensor            |
| 23 | R LL ROBOT HOME             | DI 22 | P78 | HOME/OPEN             | Sensor            |
| 24 | R PM HI VAC OPEN            | DI 23 | P79 | OPEN/UNKNOWN          | Sensor            |
| 25 | R PM HI VAC CLOSE           | DI 24 | P80 | CLOSE/UNKNOWN         | Sensor            |
| 26 | R PM ROUGHING ON            | DI 25 | P81 |                       |                   |
| 27 | SPARE                       | DI 26 | P82 |                       |                   |
| 28 | SPARE                       | DI 27 | P83 |                       |                   |
| 29 | R PM CRYO LOW TEMP          | DI 28 | P84 | LOW/UNKNOWN           | Cryo temp monitor |
| 30 | R PM CRYO HIGH TEMP         | DI 29 | P85 | HIGH/UNKNOWN          | Cryo temp monitor |
| 31 | R PM CRYO ROOM TEMP         | DI 30 | P86 | ROOM/UNKNOWN          | Cryo temp monitor |
| 32 | R PM CRYO COMPRESSOR ON     | DI 31 | P87 | ON/OFF                | Magnet relay      |

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Recipe Step Data

|                    |          |         |        |         |
|--------------------|----------|---------|--------|---------|
| ✓ Insert           | – Delete | Save    | ← Left | → Right |
|                    |          |         |        |         |
|                    |          |         |        |         |
| Step Name          | Pressure | Shutter | Mov    | Temp    |
| Time(sec)          | 10       | 30      |        | 30      |
| Pressure(100mTorr) | 5        |         |        |         |
| Ar(200sccm)        | 81       |         |        |         |
| O2(50sccm)         |          |         |        |         |
| N2(50sccm)         |          |         |        |         |
| Temperature( C )   |          |         |        | 100     |
| GUN SELECT         | -        | GUN2    | -      | -       |
| DC-3kW(PWPD301W_1  |          |         |        |         |
| DC-3kW(PWPD301W_2) |          |         |        |         |
| RF-300W(R301)      |          |         |        |         |
| Substrate Shutter  | -        | -       | -      | -       |

Step No Click 후,  
Insert Button Click

Recipe Step Data

|                    |          |      |         |         |
|--------------------|----------|------|---------|---------|
| ✓ Insert           | – Delete | Save | ← Left  | → Right |
|                    |          |      |         |         |
|                    |          |      |         |         |
| Step Name          | Pressure |      | Shutter | Mov     |
| Time(sec)          | 10       |      |         | 30      |
| Pressure(100mTorr) | 5        |      |         |         |
| Ar(200sccm)        | 81       |      |         |         |
| O2(50sccm)         |          |      |         |         |
| N2(50sccm)         |          |      |         |         |
| Temperature( C )   |          |      |         |         |
| GUN SELECT         | -        | -    | GUN2    | -       |
| DC-3kW(PWPD301W_1  |          |      |         |         |
| DC-3kW(PWPD301W_2) |          |      |         |         |
| RF-300W(R301)      |          |      |         |         |
| Substrate Shutter  | -        | -    | -       | -       |

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Recipe Step Data

|                    |          |         |        |         |
|--------------------|----------|---------|--------|---------|
| ✓ Insert           | – Delete | Save    | ← Left | → Right |
|                    |          |         |        |         |
|                    |          |         |        |         |
| Step Name          | Pressure | Shutter | Mov    | Temp    |
| Time(sec)          | 10       | 30      |        | 30      |
| Pressure(100mTorr) | 5        |         |        |         |
| Ar(200sccm)        | 81       |         |        |         |
| O2(50sccm)         |          |         |        |         |
| N2(50sccm)         |          |         |        |         |
| Temperature( C )   |          |         |        | 100     |
| GUN SELECT         | -        | GUN2    | -      | -       |
| DC-3kW(PWPD301W_1  |          |         |        | 100     |
| DC-3kW(PWPD301W_2) |          |         |        |         |
| RF-300W(R301)      |          |         |        |         |
| Substrate Shutter  | -        | -       | -      | OPEN    |

Step No Click 후  
Left or Right Button click

Recipe Step Data

|                    |          |      |          |         |
|--------------------|----------|------|----------|---------|
| ✓ Insert           | – Delete | Save | ← Left   | → Right |
|                    |          |      |          |         |
|                    |          |      |          |         |
| Step Name          | Shutter  | Mov  | Pressure | Temp    |
| Time(sec)          | 30       | 10   |          | 30      |
| Pressure(100mTorr) |          | 5    |          |         |
| Ar(200sccm)        |          | 81   |          |         |
| O2(50sccm)         |          |      |          |         |
| N2(50sccm)         |          |      |          |         |
| Temperature( C )   |          |      |          | 100     |
| GUN SELECT         | GUN2     | -    | -        | -       |
| DC-3kW(PWPD301W_1  |          |      |          |         |
| DC-3kW(PWPD301W_2) |          |      |          |         |
| RF-300W(R301)      |          |      |          |         |
| Substrate Shutter  | -        | -    | -        | -       |

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