

Gene Cloning(subcloning)

2013. 12. 31


Il Shin Kim

UNIST Central Research Facilities (UCRF)





Organization of UCRF



1. UCRF 홈페이지 접속

- UCRF homepage : <http://ucrf-eng.unist.ac.kr/main/main.php>
- 자료마당 → 자료실 → 실험기기 사용 신청서 다운로드

The screenshot shows the UCRF homepage with a navigation menu at the top. The '자료마당' (Data Market) menu is highlighted with a red box. Below the navigation, there are three main service categories: NANOFABRICATION LAB Packaging, NANOFABRICATION LAB Photolithography, and NANOFABRICATION LAB Thinfilm Deposition. A sidebar on the right contains a '자료실' (Data Room) dropdown menu with options for '자료실', 'Q & A', '관련 사이트', and '홈페이지 개선사항'. Below the main content, there are sections for '공지사항' (Notice), '자료실' (Data Room), and '교육 및 세미나' (Education and Seminar), each with a list of items and dates. A '바로가기' (Quick Links) section at the bottom right includes icons for '장비예약신청' (Equipment Reservation), '견학신청' (Field Trip), '나의장비예약현황' (My Equipment Reservation Status), '보유장비' (Owned Equipment), 'IACUC', and 'IRB'.

The screenshot shows the '자료실' (Data Room) page. At the top right, there are social media icons and a search bar. Below the navigation, there are search filters for '전체' (All) and '제목' (Title). The main content is a table of articles with columns for '번호' (No.), '분류' (Category), '제목' (Title), '작성자' (Author), '작성일' (Date), and '조회수' (Views). The table lists 62 articles in total, with the current page showing 4 out of 6 pages. Two articles are highlighted with red boxes: '기술지원 신청서 별도 양식' (Technical Support Application Form - Separate Form) and '기술지원 신청서(Technical Support) & 비용 산정 목록' (Technical Support Application Form & Cost Estimation List).

번호	분류	제목	작성자	작성일	조회수
27	기기분석실	FIB 분석 의뢰서 (Written Request for FIB Analysis)	이선이	2014-01-10	19
28	기기분석실	HR TEM 분석 의뢰서 (Written Request for HR-TEM Analysis)	박수현	2013-11-06	332
29	기기분석실	DSC, TGA, SDT, DMA, Rheometer 분석 의뢰서 (Analysis Request Form)	이경애	2013-07-19	475
30	전체	연구그룹 기자재 목록	이경선	2013-04-30	474
31	기기분석실	xps 샘플의뢰서(xps sample submit form)	ggarbi73	2012-10-25	1250
26	기기분석실	Micro-Raman self-user 교육 자료	이혜나	2013-12-18	87
25	기기분석실	Fluorometer와 spectrofluorometer self-user 교육 자료	이혜나	2013-12-18	59
24	기기분석실	FT-IR self-user 교육 자료	이혜나	2013-12-18	67
23	기기분석실	Rheometer 교육자료	이경애	2013-12-12	97
22	기기분석실	DMA 교육자료	이경애	2013-12-12	79
21	기기분석실	TGA, SDT 교육자료	이경애	2013-12-12	86
20	기기분석실	DSC 교육자료	이경애	2013-12-12	87
19	기기분석실	HPXRD 측정 및 분석방법	박지훈	2013-10-11	281
18	생체효능검증센터	실험동물 구입 신청서		2013-09-09	261
17	전체	(최종) 이용수가표	유혜정	2013-08-16	504
16	생체효능검증센터	기술지원 신청서 별도 양식	이윤진	2013-07-10	347
15	생체효능검증센터	기술지원 신청서(Technical Support) & 비용 산정 목록	이윤진	2013-07-10	342



Application of Technical Support

2. 기술 지원 신청서 작성(UNIST IVRC ANI #4)

UNIST IVRC ANI #4
No. _____

Experimental Analysis(Technical Support)

1. Applicant

Name	Office Tel./ C.P.	/
Affiliation (Department/Lab)	E-mail	
Principal Investigator (Professor in charge)	Position (Mark V)	<input type="checkbox"/> Professor <input type="checkbox"/> Researcher <input type="checkbox"/> Graduate student <input type="checkbox"/> Undergraduate student
(Sign)	Office Tel./ C.P.	/
Project Title (ACUC)	E-mail	
Date (YYYY MM DD)	IACUC Approval Number	
Species/Strains	Site	
비고	Age/ Sex/ Number	weeks/ /

2. List of Technical Support

Contents	Price (won)	Mark	항목 (원)	수기 (원)	Mark	
Inj.	Oral	1,500	Basic Experiment	Handling	500	
	Intraperitoneal	1,000		Body Weight	500	
	Subcutaneous	1,000		Tumor model	30,000/str	
	Intravenous	2,000		Other Experimental Support	30,000/str	
Blood Collection	Whole	2,000	Embryo Cryopreservation	Anatomy Training	50,000	
	Orbital venous plexus	2,000		Storage (Embryo 20days, 1 strain)	10,000/Mon	
	Biopsy(High)	5,000		Supply(basic 2 pellets)	50,000/Heads	
	Biopsy(Middle)	3,000		Embryo Freezing & Transfer	Embryo transfer	700,000
Biopsy(Low)	2,000	with Microbial monitoring (embryonic, without dipping coats)	800,000			
Anesthesia	Euthanasia	2,000	Gene Cloning Support	Embryo transfer	700,000	
	Isotranean (100m)	50,000/ea		Production of transgenic mouse vector	5,000,000	
	Zoletile	35,000/ea		Microbiological monitoring Examination (BRSB)	Virus	15,000
	Roupenin	20,000/ea			Bacteria St. Fungi	15,000
			Production of cell based expression vector	70,000		
			Production of transgenic mouse vector	5,000,000		
			Embryo transfer	700,000		
			with Microbial monitoring (embryonic, without dipping coats)	800,000		
			Embryo transfer	700,000		
			Production of transgenic mouse vector	5,000,000		
			Virus	15,000		
			Bacteria St. Fungi	15,000		
			Production of cell based expression vector	70,000		
			Production of transgenic mouse vector	5,000,000		
			Embryo transfer	700,000		
			with Microbial monitoring (embryonic, without dipping coats)	800,000		

* Embryo Freezing & Transfer - Additional Form - UNIST IVRC ANI #4-1
 ** Gene Cloning Support - Additional Form - UNIST IVRC ANI #4-2
 ☞ Reception of Inquiries: Sooh Park (sapark@unist.ac.kr, T. 052-217-5212)

Consent of providing personal information and application
 On other collect personal information such as e-mail contacts to forward the notice. After activating the purpose of the collection and use of personal information, to provide, this information will be displayed timely. You can refuse to provide personal information and this can be withdrawn at any time after consent. However, when you do not agree to provide personal information, the information provided may be limited.
 I was aware of the above and agree to use personal information.

20 Applicant : (Sign)
 I accept the above request.
 20 Representative(Receipt) : (Sign)

UNIST In Vivo Research Center.

- 의뢰자 정보기입

의뢰자의 정보와 IACUC Number 작성

- 유전자 정보 기재

유전자 Sample 상태 기재

- 기술지원 신청일 및 의뢰자 sign

3. 기술 지원 신청서 작성(UNIST IVRC ANI #4-2) 및 담당자와 조율

- 신청서 작성
- 기술지원에 관한 세부사항은 미리 담당자와 조율할것

유전자 클로닝 지원 신청서 UNIST IVRC ANI #4-2

No. _____


A Letter of Request

Approval		Client	Manager
Ref No.		Data No.	
Type of cloning	<input type="checkbox"/> Sub cloning <input type="checkbox"/> Knock-out vector cloning	<input type="checkbox"/> Mutagenesis <input type="checkbox"/> Etc ()	
University / Company			
Client	Laboratory (Supervisor)	Position	
	Name	Contact number	
	E-mail		
Gene species (Accession No.)	<input type="checkbox"/> Human <input type="checkbox"/> Rat <input type="checkbox"/> Xenopus	<input type="checkbox"/> Mouse <input type="checkbox"/> Bovine <input type="checkbox"/> Etc ()	
Kind of Vector	<input type="checkbox"/> Expression vector() <input type="checkbox"/> Knock-out vector() <input type="checkbox"/> T-Vector <input type="checkbox"/> Etc ()	No. of samples	
Sample Condition	<input type="checkbox"/> plasmid (ug/ul) <input type="checkbox"/> E.coli <input type="checkbox"/> Etc ()		
Misc. Information			

UNIST In Vivo Research Center. IVRC

- 서명
- Type of Cloning 기재
- 기술지원 신청자의 정보를 기재
- 유전자 정보 기재
- 유전자 Sample 상태 기재

Gene Cloning의 기본적 흐름




Defining Cloning

◆ “Cloning” is a loaded term that can be used to mean very different things.

- Cutting a piece of DNA from one organism and inserting it into a vector where it can be replicated by a host organism. (Sometimes called subcloning, because only part of the organism's DNA is being cloned.)

- Using nuclear DNA from one organism to create a second organism with the same nuclear DNA





DNA Cloning

DNA cloning is a technique for reproducing DNA fragments. It can be achieved by two different approaches:

- cell based
- using polymerase chain reaction (PCR).

a vector is required to carry the DNA fragment of interest into the host cell.

DNA cloning allows a copy of any specific part of a DNA (or RNA) sequence to be selected among many others and produced in an unlimited amount.

This technique is the first stage of most of the genetic engineering experiments:

- production of DNA libraries
- PCR
- DNA sequencing


Massive amplification of DNA sequences

Stable propagation of DNA sequences

A single DNA molecule can be amplified allowing it to be:

- Studied - Sequenced
- Manipulated - Mutagenised or Engineered
- Expressed - Generation of Protein


Cloning a Piece of DNA



Cut plasmid vector
with EcoR I



Excise DNA insert of interest from source
using EcoR I



Ligate the insert of interest
into the cut plasmid



Plasmid Cloning Strategy

◆ Involves five steps:

1. Enzyme restriction digest of DNA sample.
2. Enzyme restriction digest of DNA plasmid vector.
3. Ligation of DNA sample products and plasmid vector.
4. Transformation with the ligation products.
5. Growth on agar plates with selection for antibiotic resistance.


※ Strategy depends on the starting information and desired endpoint.

Starting Information or Resources:


- Protein sequence
- Positional cloning information
- mRNA species / sequence
- cDNA libraries
- DNA sequence known or unknown
- Genomic DNA libraries
- PCR product

Plasmid Cloning Strategy


Step 1. Re Digestion of DNA Sample



Step 2. Re Digestion of Plasmid DNA



Step 3. Ligation of DNA Sample and Plasmid DNA




Step 4. Transformation of Ligation Products


The process of transferring exogenous DNA into cells is called “transformation”

There are basically two general methods for transforming bacteria. The first is a chemical method utilizing CaCl_2 and heat shock to promote DNA entry into cells.

A second method is called electroporation based on a short pulse of electric charge to facilitate DNA uptake.





Step 5. Growth on Agar Plates





Plasmid Vector

- Plasmid vectors are $\approx 1.2\text{--}3\text{kb}$ and contain:
 - replication origin (ORI) sequence
 - a gene that permits selection,
 - Here the selective gene is *amp^r*; it encodes the enzyme b-lactamase, which inactivates ampicillin.
 - Exogenous DNA can be inserted into the bracketed region .
- Many cloning vectors contain a **multiple cloning site** or **polylinker**: a DNA segment with several unique sites for restriction endo- nucleases located next to each other
- Restriction sites of the polylinker are not present anywhere else in the plasmid.
- Cutting plasmids with one of the restriction enzymes that recognize a site in the polylinker does not disrupt any of the essential features of the vector



Plasmid Vectors

- ◆ Plasmids are circular pieces of DNA found naturally in bacteria.
- ◆ Plasmids can carry antibiotic resistance genes, genes for receptors, toxins or other proteins.
- ◆ Plasmids replicate separately from the genome of the organism.
- ◆ Plasmids can be engineered to be useful cloning vectors.
- ◆ Plasmid vectors can be designed with a variety of features:
 - Antibiotic resistance
 - Colorimetric “markers”
 - Strong or weak promoters for driving expression of a protein





The cloning marker for this plasmid is the *lacZ* gene.

Restriction Enzymes

Restriction Enzymes (also called Restriction Endonucleases) are proteins that cleave DNA molecules at specific sites, producing discrete fragments of DNA.

Generally restriction enzymes are thought to protect bacterial cells from phage (bacterial virus) infection. Bacterial cells that contain restriction enzymes can “cut up” invasive viral DNA without damaging their own DNA.

- Restriction enzymes: enzymes that cut DNA in specific places function:
- Inactivate foreign DNA
- Breaks only palindrome sequences, i.e. those exhibiting two-fold symmetry
- Important in DNA research, i.e. sequencing, hybridization
- Companies purify and market restriction enzymes






Isolating Genes

- ◆ Herbert Boyer and Stanley Cohen built on the work of Berg, Nathans and Smith to use restriction enzymes to isolate a single gene, place it into a plasmid vector.
- ◆ Bacterial cells were then transformed with the recombinant plasmid.
- ◆ The bacteria host cells replicated the plasmid, producing many copies of the gene, thus amplifying it.
- ◆ The practical application was that expensive human protein products, like insulin, which were used to treat disease, could eventually be produced from recombinant molecules in the laboratory using bacteria or another host.
- ◆ Human protein products like insulin could be used in very large quantities from the recombinant molecule. Patients no longer had to use insulin isolated from pigs or cows.



Performing the Restriction Digests

- ◆ You will need to set up a restriction digest of your plasmid vector and your DNA of interest
- ◆ Restriction enzymes all have specific conditions under which they work best. Some of the conditions that must be considered when performing restriction digest are: temperature, salt concentration, and the purity of the DNA



Purify your DNA Fragments


- ◆ The insert of interest that you want to clone into your plasmid needs to be separated from the other DNA
- ◆ You can separate your fragment using Gel Electrophoresis
- ◆ You can purify the DNA from the gel by cutting the band out of the gel and then using a variety of techniques to separate the DNA from the gel matrix



Gel Elution

Ligation

- ◆ Ligation is the process of joining two pieces of DNA from different sources together through the formation of a covalent bond.
- ◆ DNA ligase is the enzyme used to catalyze this reaction.
- ◆ DNA ligation requires ATP.





Transformation in Bacteria cells

1. Transforming Bacteria

- ◆ After you create your new plasmid construct that contains your insert of interest , you will need to insert it into a bacterial host cell so that it can be replicated.
- ◆ The process of introducing the foreign DNA into the bacterial cell is called transformation.

2. Competent Host Cells

- ◆ Not every bacterial cell is able to take up plasmid DNA.
- ◆ Bacterial cells that can take up DNA from the environment are said to be competent.
- ◆ Can treat cells (electrical current/divalent cations) to increase the likelihood that DNA will be taken up
- ◆ Two methods for transforming: heat shock and electroporation

3. Selecting for Transformants

- ◆ The transformed bacteria cells are grown on selective media (containing antibiotic) to select for cells that took up plasmid.
- ◆ For blue/white selection to determine if the plasmid contains an insert, the transformants are grown on plates containing X-Gal and IPTG.



Expressing your cloned gene

- ◆ Even if your plasmid contains insert, it may not be able to generate functional protein from your cloned DNA.
 - The gene may not be intact, or mutations could have been introduced that disrupt it.
 - The protein encoded by the gene may require post-translational modifications (i.e., glycosylation or cleavage) to function.
 - Also, some enzymes are a complex of peptides expressed from separate genes.

- ◆ Expression of a cloned gene can be accomplished by:
 - The *E. coli* host
 - Mammalian cells (if the plasmid used is designed for expression in mammalian cells)
 - Using an in vitro using a cell-free system.

Cloning into a Plasmid

