
UNIST
Ulsan National Institute of
Science and Technology

Introduction of LABs

Be one of us! Be the best of all!

Contents

Introduction	04
A message from the president	05
Organization	04
A List of Instruments	05
Electron Microscopy	06
X-Ray Analysis	12
Spectroscopic Analysis	16
Surface Analysis	22
Mass Spectrometry	29
Thermal Analysis	35
Environmental Analysis	38
A List of Instruments (Nanofabrication Lab)	41
Photolithography	42
Etching Process	46
Thin Film Deposition	50
Measurement	58
Packaging	62
A List of Instruments (Machine Shop)	41
CNC Machine Tool	64
General purpose Machine Tool	68
Measuring Instrument	72
Procedure to use	78
How to get to UNIST	79

A Message from the President

UNIST has been founded with the aim of becoming not only one of Korea's most prominent national universities but also one of leading global universities in education and research of science and technology. UNIST is located in Ulsan Metropolitan City where Korea's leading industries including automobile, ship building and petroleum chemistry are clustered. With this location, UNIST has the potential of emerging as a hub of industry-academic cooperation in high technology. One of UNIST's long term visions is being a leading global university by focusing on interdisciplinary studies of science and technology. By recruiting highly qualified faculty, staff and students from all over the world, UNIST pledges to provide an optimal educational environment to our students. research activities will be fo-

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cused on some interdisciplinary areas which can generate basic principles and on technological innovation. And, coursework for graduate students will mainly consist of various research activities. Construction of a first-class ubiquitous campus and other infrastructure was opened in 2009. We plan to e promote education and research activities. UNIST is aimed at serving as a base for cultivating global leaders who can lead our society, industry and the academic world, and creating innovative theories and cutting-edge technologies that can make a great impact on our future society. As a first president, I will be committed to develop UNIST as a world's leading university in the 21st century. UNIST has been founded with the aim of becoming not only one of Korea's most prominent national universities but also one of the leading global universities in education and research of science and technology.

UNIST is located in Ulsan Metropolitan City where Korea's leading industries including automobile, ship building and petroleum chemistry are clustered. With this location, UN-

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IST has the potential of emerging as a hub of industry-academic cooperation in high technological innovation. One of UNIST's long-term visions is being a leading global university by focusing on interdisciplinary studies of in science and technology. By recruiting highly qualified faculty, staff and students from all over the world, UNIST pledges to provide an optimal educational environment to our students. Research activities will be focused on some interdisciplinary areas which can generate basic principles and on technological innovation. And, coursework for graduate students will mainly consist of various research activities. Construction of Our first-class ubiquitous campus and other infrastructure was opened in 2009. We plan to e-promote education and research activities. UNIST is aimed at serving as a base for cultivating global leaders who can lead our society, industry and the academic world, and creating innovative theories and cutting-edge technologies that can make a great impact on our future society. As a the first president, I will be committed to develop UNIST as a world's leading university in the 21st century.

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

Facility

Lab	Total	Nanofabrication	Machine Shop	Environmental Analysis	Biomed Imaging	Analysis
Area (m ²)	5,940	730	1,700	230	200 (공학 1동)	3,080

A List of Instruments

Analysis Lab

Lab	Equipment	Model	Maker	Page
Electron Microscopy	High Resolution Transmission Electron Microscope	JEM-2100F (Cscorrector)	JEOL, Japan	6
	Transmission Electron Microscope	JEM-2100	JEOL, Japan	7
		JEM-1400	JEOL, Japan	
	Field Emission Scanning Electron Microscope	Quanta 200	FEI, USA	8
		Nano230	FEI, USA	9
	Dual-beam Focused Ion Beam	Quanta 3D FEG	FEI, USA	10
X-Ray Analysis	High Power X-Ray Diffractometer	D/Max2500V/PC	Rigaku, Japan	11
	High Resolution X-Ray Diffractometer	D8 Discover	Bruker, Germany	12
	X-Ray Fluorescence Spectrometer	T8 Tiger	Bruker, Germany	13
	Single Crystal X-Ray Diffractometer	R-AXIS RAPID II	Bruker, Germany	14
Spectroscopic Analysis	FT-NMR (600 MHz)	VNMRS 600	Varian, USA	15
	ICP-OES	720-ES	Varian, USA	16
	Fluorometer	Cary Eclipse	Varian, USA	17
	UV-Vis-NIR	Cary 5000	Varian, USA	18
		670	Varian, USA	19
		620	Varian, USA	
	Sub-micron size & Zeta potential measuring system	Nano ZS	Malvern, UK	20
Surface Analysis	X-ray Photoelectron Spectroscopy	K-alpha	Thermo Fisher, UK	21
	Variable Temperature Scanning Tunneling Microscope	UHV VT 7000	RHK Technology	22
	Combined AFM-Raman Microscope	Alpha300S	WITec, USA	23
	Confocal Raman Microscope	Alpha300R	WITec, U	
	Atomic Force Microscope	Multimode V	Veeco, USA	24
	Atomic Force Microscope	Dimensim 3100	Veeco, USA	
	Electrochemical Atomic Force Microscope	Agilent 5500	Agilent, USA	25
	Confocal Laser Scanning Microscope	OLS3100	Olympus, Japan	26
	Cryogenic Probe Station	CRX-4K	Lakeshore, USA	27
		TOF/SIMS	TOF SIMS 5	ION TOF, Germany
Mass Analysis	MALDI TOF/TOF	Ultraflex III	Bruker, Germany	29
	GC/MS/MS	450-GC	Varian, USA	30
		320-MS	Varian, USA	
	LC/MS/MS	HCT Basic System	Bruker, Germany	31
	Gel permeation chromatography (GPC)/ Multi Angle Light Scattering (MALS)	Agilent 1200 series	Agilent, USA	32
		miniDAWN TREOS	Wyatt, USA	
Element Analyzer	Flash 2000	Thermo, Netherlands	33	
Thermal Analysis	Thermogravimetric Analyzer	Q400	TA Instruments, USA	34
	Differential Scanning Calorimetry	Q600	TA Instruments, USA	35
	Dynamic Mechanical Analysis	Q800	TA Instruments, USA	36
Environmental Analysis	HRGC/HRMS	Autospec Premier	Waters, USA	37
	ICPMS	ELAN DRC-e	Perkin Elmer, USA	38
		7890A GC	Agilent, USA	39
	GC/MSD	5975C MSD	Agilent, USA	
	GC/ECD	7890A GC	Agilent, USA	40

High Resolution Transmission Electron Microscope

Model

- JEM-2100F (Cs corrector) (JEOL, JAPAN)

Specifications

- Acceleration voltage : 200 kV
- Lattice resolution : 0.102 nm
- BF-STEM resolution : 0.14 nm
- HAADF-STEM resolution : 0.096 nm
- Probe stability : 1 nm/min or less
- Stage stability : 1 nm/min or less
- SEDS resolution : 128 eV
- EELS resolution : 0.8 eV
- Contamination : 1 nm/min or less

Manager : Kim Young-Ki

- e-mail : cclock95@unist.ac.kr
- Phone : 052-217-4033

Transmission Electron Microscope

Model

- JEM-2100 (JEOL, JAPAN)
- JEM-1400 (JEOL, JAPAN)

Specifications

- Acceleration voltage :
JEM-2100 : 200 kV / JEM-1400 : 120 kV
- Lattice resolution :
JEM-2100 : 0.23 nm / JEM-1400 : 0.38 nm
- Point resolution :
JEM-2100 : 0.14 nm / JEM-1400 : 0.20 nm
- Specimen tilting :
- JEM-2100 : X = $\pm 35^\circ$, Y = $\pm 30^\circ$
- JEM-1400 : X = $\pm 25^\circ$, Y = $\pm 70^\circ$
- Image recording system : CCD
- EDS resolution : 132 eV

Manager : Kim Young-Ki

- e-mail : cclock95@unist.ac.kr
- Phone : 052-217-4033

Basic Principle

HR-TEM is equipped with field emission electron gun featuring high brightness and long-term stability. Nano scale observation of ultra microstructures at atomic level and high sensitivity analysis of ultra microstructures is possible. It is possible to get atomic resolution images as well as spectrum images in the nano region using EELS.

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Basic Principle

TEM is used to analyze micro-structures of materials with high spatial resolution. The high voltage electron beam generated from electron gun is illuminated on thin film specimen. The beam penetrating the specimen passes through an array of magnetic lenses and forms a high resolution electron image of the electron diffraction pattern.

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Field Emission-Scanning Electron Microscope

Model

- Quanta 200 (FEI, USA)

Specifications

- Schottky type thermal FE gun
- Resolution :
 - < 1.2 nm @ 30 kV high vacuum mode
 - < 1.5 nm @ 30 kV low vacuum mode
 - < 1.5 nm @ 30 kV ESEM mode
- Beam current : > 100 nA
- Chamber vacuum : Up to 40 mbar
- 4-axis motorized stage
- Attachment : EDS
- 1000 or 1500 °C heating stage
- Cooling stage (-20 °C)

Manager : Lim Dong-Ju

- e-mail : djlim@unist.ac.kr
- Phone : 052-217-4032

Field Emission-Scanning Electron Microscope

Model

- NanoSem 230 (FEI, USA)

Specifications

- Schottky type thermal FE gun
- Conventional operation :
 - < 1.0 nm @ 15 kV
 - < 1.6 nm @ 1 kV
- Reduced charge & contamination operation
- Probe current : 0.6 pA ~ 100 nA
- 4-axis motorized stage
- Beam deceleration standard
- Attachment : EDS

Manager : Lim Dong-Ju

- e-mail : djlim@unist.ac.kr
- Phone : 052-217-4032

Basic Principle

SEM is a microscope that uses electrons instead of light to form an image. The scanning electron microscope has many advantages over traditional microscopes. The SEM has a large depth of field, which allows more of a specimen to be in focus at one time. The SEM also has much higher resolution, so closely spaced specimens can be magnified at much higher levels. Because the SEM uses electromagnets rather than lenses, the researcher has much more control in the degree of magnification. Field-emission scanning electron microscope provides ultra high resolution image down to 1nm resolution thanks to the inherent brightness of field – emission electron gun.

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Basic Principle

A finely focused electron beam scanned across the surface of the sample generates secondary electrons, backscattered electrons, and characteristic X-rays. These signals are collected by detectors to form images of the sample displayed on a cathode ray tube screen. Elemental composition of features seen in the SEM image can then be immediately analyzed by EDS.

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- e-mail : djlim@unist.ac.kr
- Phone : 052-217-4032

Basic Principle

The Dual Beam FIB has a multitude of capabilities including high resolution electron imaging, ion imaging, nano device fabrication and material deposition. Simultaneous patterning and imaging are possible to nm resolution for both imaging and machining. The FIB with both ion and electron beams enables high-precision cross-section processing and TEM/3DAP sample preparation.

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- e-mail : sophia@unist.ac.kr
- Phone : 052-217-4035

Basic Principle

X-ray diffraction is a non-destructive analytical technique which reveals information about the crystallographic structure, chemical composition, and physical properties of materials and thin films. X-ray diffraction is based on constructive interference of monochromatic x-rays and a crystalline sample. The High Power XRD is equipped with a high-temp attachment for temperature variation and a battery cell for charging and discharging.

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High Resolution X-Ray Diffractometer

Model

- D8 Advance (Bruker, Germany)

Specifications

- X-ray source : rotating anode
- Focus : line focus
- X-ray generator : single, 220 V, 3 kW
- θ - 2θ based goniometer
- Angular range (θ) : 0 ~ 360°
(2θ) : - 110 ~ 168° more

Manager : Park Ji-Hye

- e-mail : sophia@unist.ac.kr
- Phone : 052-217-4035

Basic Principle

X-ray diffraction is a non-destructive analytical technique which reveals information about the crystallographic structure, chemical composition and physical properties of materials and thin films. It is based on observation of the scattered x-ray beam after hitting a sample as a function of incident and scattered angle, polarization and wavelength or energy.

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- e-mail : sophia@unist.ac.kr
- Phone : 052-217-4035


Basic Principle

X-ray fluorescence is the emission of characteristic "secondary" x-rays from a material that has been excited by bombarding with high-energy x-rays or gamma rays. XRF is able to analyze quality and quantity of various elements of solid materials such as geological specimens, slugs, soils, ceramics very quickly as well as simultaneously. Qualitative analysis of elements ranging from 4Be to 92U is made and quantitative analysis of 10 major elements in the case of geological sample is usually made in concentration of 100 ppm – 100 %. With the advantage of non-destructive analysis, precision and accuracy by repetitive measurements is certified. Moreover, sample preparation is simple.

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Application

- Element analysis 원소분석
- Semi quantitative analysis 반정량분석

Single Crystal X-Ray Diffractometer

Model

- R-AXIS RAPID II (Rigaku, Japan)

Specifications

- X-ray generator with high stability :
3kw, max : 60 KV/ 60 mA
- X-ray tube : Mo source
- Goniometer : decoupled theta and 2 theta
- Analyzing crystal
- Optics and chamber
- Detectors
- Imaging Plate(IP) detector : + 5 ~ + 35 °C

Manager : Kim Young-Ki

- e-mail : cclock95@unist.ac.kr
- Phone : 052-217-4033

Basic Principle

R-AXIS RAPID II is 2 demension XRD which has curved Imaging plate. It can measure rapidly a variety of samples such as small parts of bulk, normal powder, film and fiber samples. It also can distinguish single crystal from multi crystal indicating the difference of particle size and alingment through the 2D image patterns. 2D image X_ray Patterns from 2D detector can tell us the information of samples' alignments or crystal structures. Parallel beam optical system is unconfined by the form of sample and it can offer various sample holders such as powder, bulk, film, fiber, trace sample, liquid and thin film.

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- VNMRS 600 (Varian, USA)

Specifications

- Premium shielded 600 MHz magnet
- 3 channel console
- Available probes
 - Dual broad band auto X probe
 - Automated triple resonance probe
 - Nano TM probe
 - 5mm double resonance MAS solid probe
- Possibility of advanced techniques
 - 1D/2D/3D
 - 1H, 13C, 19F, 31P
 - COSY, HSQC, HMBC, NOESY, TOCSY

Manager : Cho Mi-Sun

- e-mail : shail019@unist.ac.kr
- Phone : 052-217-4034

Basic Principle

NMR is a technique for determining the structure of most organic and some inorganic molecules. NMR studies a magnetic nucleus by aligning it with a very powerful external magnetic field and disturbing this alignment using an electromagnetic field. Different atoms within a molecule resonate at different frequencies at a given field strength. The observation of the resonance frequencies of a molecule allows us to discover structural information of the molecule.

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- 720-ES (Varian, USA)

Specifications

- Wavelength range : 167 ~ 785 nm
- Simultaneous analysis
- 40 MHz free running RF generator
- CCD detector with image map tech.
- Cooled cone interface mechanism
- Plasma : axially viewed system

Manager : Cho Mi-Sun

- e-mail : shail019@unist.ac.kr
- Phone : 052-217-4034

Fluorometer

Model

- Cary Eclipse (Varian, USA)

Specifications

- Wavelength range
 - Excitation 200 ~ 1,100 nm
 - Emission 200 ~ 1,100 nm
- Wavelength accuracy : < 1.5 nm
- Limiting resolution : 1.5 nm
- Wavelength reproducibility : < 0.2 nm
- Source : Xe pulsed lamp
- Detector : photomultiplier tube

Manager : Cho Mi-Sun

- e-mail : shail019@unist.ac.kr
- Phone : 052-217-4034

Basic Principle

ICP-OES is a technique for element analysis of any substances in the ppb ~ ppm concentration range. It is based on the introduction of the liquid sample to be analyzed towards an argon plasma, created by a radiofrequency generator. The plasma energy is transferred to the atoms of every sample elements, promoting the excitation of their electrons to higher energy levels. These excited electrons are decayed to low energy levels and electromagnetic radiation in the UV/Vis range is emitted.

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Basic Principle

The fluorometer is capable of measuring excitation and emission spectra in the UV, Vis and NIR regions, roughly from 200 nm to 1,100 nm. Absorption of UV radiation by a molecule excites it from the electronic ground state to the excited state. Fluorescence occurs when the molecule returns to the electronic ground state from the excited state.

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- e-mail : shail019@unist.ac.kr
- Phone : 052-217-4034

Basic Principle

UV-Vis-NIR is used for optical absorbance, transmittance and reflectance measurements in the wavelength range (175 ~ 3,300 nm). In principle the technique is similar to IR-absorbance. When a sample of an unknown compound is exposed to light, certain functional groups within the molecule absorb light of different wavelengths in the UV or Vis or NIR region. It is used for qualitative and quantitative analysis of materials.

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Basic Principle

The FT-IR measures the IR absorbance or reflectance spectrum to determine vibrational frequencies, which provide information on chemical structure and bonding. The spatial resolution of this technique is in the order of 10 μm , enabling it to determine the composition of very small features on various substrates. The technique can also study dopants in glasses and semiconductors.

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Sub micron size & Zeta potential measuring system

Model

- Nano ZS (Malvern, UK)

Specifications

- Power : single 220 V, 60 Hz
- Particle measurement
 - volume range of 12 μl ~ 1.0 mL
 - sample concentration > 40 % by weight
 - min. concentration 0.1 mg/mL
- Zeta potential measurement
 - Principle : laser doppler velocimetry
 - Range : - 150 mV ~ + 150 mV

Manager : Park Ji-Hye

- e-mail : sophia@unist.ac.kr
- Phone : 052-217-4035

X - ray Photoelectron Spectroscopy

Model

- K-alpha (Thermo Fisher, UK)

Specifications

- Double-focusing hemispherical analyzer
- Multi-element, high-transmission spectrometer input lens
- 128-channel detector for high quality
- Energy range 200 eV to 3 keV
- High flux even at low beam energy
- 250 mm Rowland circle monochromator

Manager : Lim Dong-Ju

- e-mail : djlim@unist.ac.kr
- Phone : 052-217-4032

Basic Principle

Particles, emulsions and molecules in suspension undergo Brownian motion. If the particles or molecules are illuminated with a laser, the intensity of the scattered light fluctuates at a rate that is dependent upon the size of the particles as smaller particles are kicked further by the solvent molecules and move more rapidly. Analysis of these intensity fluctuations yields the velocity of the Brownian motion and hence the particle size using the Stokes-Einstein relationship. Dynamic light scattering is a non-invasive and well-established technique for measuring the size of molecules and particles typically in the submicron region.

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Basic Principle

X-ray photoelectron spectroscopy (XPS) is a quantitative spectroscopic technique that measures the elemental composition, empirical formula, chemical state and electronic state of the elements that exist within a material. XPS spectra are obtained by irradiating a material with a beam of aluminum or magnesium X-rays while simultaneously measuring the kinetic energy (KE) and number of electrons that escape from the top 1 to 10 nm of the material being analyzed.

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Variable Temperature Scanning Tunneling Microscope

Model

- UHV VT 7000 (RHK Technology, USA)

Specifications

- Scan range (X, Y) : 5 μm
- Resolution (X, Y) : 0.5 \AA
- Resolution (Z) : 0.05 \AA
- Thermal drift (X, Y & Z) < 1 $\text{\AA}/\text{min}$
- Cooling temperature : 30 K or lower
- Heating temperature : 1000 $^{\circ}\text{C}$ or higher

Manager : Lim Dong-Ju

- e-mail : djlim@unist.ac.kr
- Phone : 052-217-4032

Basic Principle

STM probes the density of states of a material using tunneling current. For STM, good resolution is considered to be 0.05 nm lateral resolution and 0.005 nm depth resolution. The STM can be used not only in ultrahigh vacuum but also in air and various other liquid or ambient gas and at temperatures ranging from near zero Kelvin to a thousand degrees Celsius. The STM is based on the concept of quantum tunneling. When a conducting tip is brought very near to a metallic or semiconducting surface, a bias between the two can allow electrons to tunnel through the vacuum between them. Variations in current as the probe passes over the surface are translated into an image. STM can be a challenging technique, as it requires extremely clean surfaces and sharp tips.

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Combined AFM & Confocal Raman Microscope

Model

- Alpha300S (WITec, Germany)
- Alpha300R (WITec, Germany)

Specifications

- Scanning near-field optical microscope
- AFM mode
 - Acoustic AC mode
 - Contact mode
- Laser
 - Extension confocal Raman to 785 nm
 - Frequency doubled Nd : YAG laser
 - Raman laser coupler 633 nm
 - Helium Neon (HeNe) laser : 532 nm
- Liquid imaging extension
- Heating and freezing cryostat

Manager : Park Ji-Hye

- e-mail : sophia@unist.ac.kr
- Phone : 052-217-4035

Basic Principle

By simply rotating the objective turret, the user can choose a mode among confocal microscopy, SNOM or AFM. This allows topography and optical images to be acquired simultaneously. Also it is possible to start with confocal Raman microscopy and combine to atomic force microscopy or vice versa. With such a combination, chemical information can be directly linked to structural AFM information from the same sample area using only one instrument.

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- Multimode V (Veeco, USA)
- Dimensim 3100 (Veeco, USA)

Specifications

- Resolution
 - Noise level RMS : < 0.3 Å RMS
 - Vertical resolution : 0.1 Å
 - Lateral resolution : 1 Å
- Controller
 - X,Y,Z axis : 16 bit, Total : 48 bit
 - High-speed DSP/ADC
 - 2 high-speed ADCs (50 MHz 14 bit)

Manager : Lim Dong-Ju

- e-mail : djlim@unist.ac.kr
- Phone : 052-217-4032

Electrochemical Atomic Force Microscope

Model

- Agilent 5500 (Agilent, USA)

Specifications

- 100 μm x 100 μm scanning range with precision, closed-loop nano positioning
- 0.05 % non-linearity
- 0.5 nm RMS X-Y position noise
- Capacitive closed-loop sensors
- Compatible with all imaging modes, multipurpose scanners & liquid cells
- PicoLith software

Manager : Lim Dong-Ju

- e-mail : djlim@unist.ac.kr
- Phone : 052-217-4032

Basic Principle

AFM or SPM is a branch of microscopy that forms images of surfaces using a physical probe that scans the specimen. An image of the surface is obtained by mechanically moving the probe in a raster scan of the specimen line by line and recording the probe-surface interaction as a function of position.

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Basic Principle

Electrochemical SPM offers a low-noise potentiostat / galvanostat for in situ EC-SPM studies – both EC-STM and EC-AFM. When combined with temperature control, it is possible to obtain valuable information about electrochemical processes that would otherwise be inaccessible. The addition of environmental control allows imaging with no dissolved oxygen in either aqueous or non-aqueous solutions.

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Confocal Laser Scanning Microscope

Model

- OLS3100 (Olympus, Japan)

Specifications

- Power : single 220 V, 4 kW
- Operation temperature : 22°C ± 3°C
- Microscopy : conventional/laser mode
- Detector
 - 2 PMT fluorescence detectors
 - 2 individual channels simultaneously

Manager : Park Ji-Hye

- e-mail : sophia@unist.ac.kr
- Phone : 052-217-4035

Basic Principle

In a confocal laser scanning microscope, a laser beam passes through a light source aperture and then is focused by an objective lens into a small focal volume within a fluorescent specimen. A mixture of emitted fluorescent light as well as reflected laser light from the illuminated spot is then collected by the objective lens. A beam splitter separates the light mixture by allowing only the laser light to pass through and reflecting the fluorescent light into the detection apparatus. After passing a pinhole, the fluorescent light is detected by a photo-detector (a photomultiplier tube or avalanche photodiode).

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- CRX-4K (Lake Shore, USA)

Specifications

- Temperature : 1.5 K ~ 475 K
- Six micro-manipulated probe arms
- 4-inch wafer probe capabilities
- Cryogen-free CCR (Closed Cycle Refrigerator)
- Vertical or horizontal field magnets
- High vacuum to 10⁻⁷ Torr
- Load-lock

Manager : Park Ji-Hye

- e-mail : sophia@unist.ac.kr
- Phone : 052-217-4035

Basic Principle

A probe station is used to physically acquire signals from the internal nodes of a semiconductor device. If the device is being electrically stimulated, the signal is acquired by the mechanical probe and is displayed on an oscilloscope. This is a versatile micro-manipulated probe station used for non-destructive testing of devices on full and partial wafers up to 2 inch in diameter. This is a platform for measurement of electrical, electro-optical, parametric, high Z, DC, RF, and microwave properties of material and test devices.

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- TOF SIMS 5 (ION TOF, Germany)

Specifications

- Mass resolution (@ 29 u) : > 10,000 (for Bi+)
- Sensitivity (@ 29 u) : > 3×10^{-8} Al⁺/nC
- Mass range : > 9,000 amu
- Base Pressure : < 5.0×10^{-10} torr
- A self-adjusting charge compensation system
- Primary source : Pulsed Bi cluster ion source
- O₂ and Cs dual sources for depth profiling
- A motorized five-axes UHV sample stage
- Two CCD high resolution video cameras

Manager : Lee Kyung-Sun

- e-mail : kslee@unist.ac.kr
- Phone : 052-217-4036

Basic Principle

Secondary ion mass spectrometry (SIMS) is a technique to analyze the composition of solid surfaces and thin films by sputtering the surface of the specimen with a focused primary ion beam and collecting and analyzing ejected secondary ions. Time-of-flight (TOF) is a method of mass spectrometry in which ions are accelerated by an electric field of known strength. The velocity of the ion depends on the mass-to-charge ratio. The time that it subsequently takes for the particle to reach a detector at a known distance is measured. This time will depend on the mass-to-charge ratio of the particle (heavier particles reach lower speeds).

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- Ultraflex III (Bruker, Germany)

Specifications

- High mass resolution at every mass from 700 – 6,000 Da
- Resolution $\geq 1,100$ & mass accuracy of ≤ 80 ppm in linear mode for protein
- Mass accuracy of ≤ 3 ppm for peptide mass range in reflectron mode
- Sensitivity in low attomole range with high S/N ratio in TOF/TOF mode

Manager : Lee Kyung-Sun

- e-mail : kslee@unist.ac.kr
- Phone : 052-217-4036

Basic Principle

Matrix-assisted laser desorption/ionization (MALDI) is a soft ionization technique used in mass spectrometry, allowing the analysis of biomolecules and large organic molecules which tend to be fragile and fragment when ionized by more conventional ionization methods. TOF/TOF is a tandem mass spectrometry method where two time-of-flight mass spectrometers are used consecutively. The first TOF-MS is used to separate the precursor ions, and the second TOF-MS analyzes the product ions. An ion gate for selecting the precursor ion, an ion fragmentation region (e.g. a collision cell) and an ion accelerator may be provided between the first and second TOF-MS.

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- 450-GC & 320-MS (Varian, USA)

Specifications

- Ionization mode
 - Electron ionization
 - Positive/negative chemical ionization
- Mass range : 10 ~ 2,000 amu
- Mass axis stability : +/-0.1 u over 24 hrs
- Resolution : 0.7 amu at 1,250 amu/sec
0.6 amu at 500 amu/sec
- Mass type : tandem triple quadrupole
- Library : NIST library

Manager : Lee Kyung-Sun

- e-mail : kslee@unist.ac.kr
- Phone : 052-217-4036

LC / MS / MS

Model

- HCT Basic System (Bruker, Germany)

Specifications

- Ion source : ESI, APCI
- Mass range : 50 ~ 3,000 (m/z)
- Resolution : 0.6 u
- Accuracy : 0.1 ~ 0.3 % absolute
- Scan speed : 26,000 amu/sec
- Mass accuracy : +/-0.15 u
- Electro sensitivity
 - MS : reserpine 5 pg S/N 10:1
 - MS/MS : reserpine 1 pg S/N 50:1

Manager : Lee Kyung-Sun

- e-mail : kslee@unist.ac.kr
- Phone : 052-217-4036

Basic Principle

GC/MS/MS makes an effective combination for chemical analysis. Mixtures of volatile substances are separated by gas chromatography. Via a transfer line the fractions of each peak are transferred into the ion source where the ionization takes place. After acceleration, the ions pass through a magnetic-field which causes deviation of their linear flight curve according to their respective mass. After this magnetic field the detector counts the ions.

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Basic Principle

LC/MS/MS is a powerful and widely used tool in qualitative and quantitative residue and contaminant analysis. The mixture is separated with an HPLC, preferably using a reversed-phase column. The analytes can be ionized with a suitable ion source by various methods, followed by partial fragmentation. After acceleration they are deflected by a magnetic field which resolves them according to their mass.

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- e-mail : shail019@unist.ac.kr
- Phone : 052-217-4034

Basic Principle

Gel permeation chromatography (GPC) is a term used for when the separation technique size exclusion chromatography (SEC). The smaller analytes can enter the pores more easily and therefore spend more time in these pores, increasing their retention time. Conversely, larger analytes spend little if any time in the pores and are eluted quickly. GPC is often used to determine the relative molecular weight of polymer. Multiangle light scattering (MALS) is a technique for determining, independently, the absolute molar mass and the average size of particles in solution, by detecting how they scatter light.

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Element Analyzer

Model

- Flash 2000 (Thermo Scientific, Netherlands)

Specifications

- Measuring elements : C, H, N, S & O
- Measuring range : 100 ppm ~ 100 %
- Sample size : 0.01 ~ 500 mg
- Accuracy : 0.1 ~ 0.3 % absolute
- Detector : TCD
- Furnace temperature : 1,100 °C
- Analysis time
 - CHNS : 10 min
 - Oxygen : 5 min

Manager : Cho Mi-Sun

- e-mail : shail019@unist.ac.kr
- Phone : 052-217-4034

Basic Principle

Element analyzer determines the percentage weights of carbon, hydrogen, nitrogen, sulfur and oxygen of a sample. This information is important to help determine the structure of an unknown compound, as well as to help prove the structure and purity of a synthesized compound. It is based on the dynamic flash combustion of the sample followed by reduction, trapping, complete GC separation and detection of the products by thermal conductivity detector.

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- e-mail : cclock95@unist.ac.kr
- Phone : 052-217-4033

Basic Principle

Thermogravimetric Analyzer (TGA) measures the amount of weight change of a material, either as a function of increasing temperature or as a function of time in an atmosphere of nitrogen, argon and air.

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Differential Scanning Calorimetry

Model

- Q600 (TAI, USA)

Specifications

- Temperature range : 725 $^{\circ}$ C
- Cooling accessories : - 180 ~ 725 $^{\circ}$ C
- Temperature accuracy : \pm 0.1 $^{\circ}$ C
- Calorimetric reproducibility : \pm 0.1 %
- Baseline curvature : 10 μ W
- Sensitivity : 0.2 μ W

Manager : Kim Young-Ki

- e-mail : cclock95@unist.ac.kr
- Phone : 052-217-4033

Basic Principle

Differential scanning calorimetry (DSC) is a thermo analytical technique in which the difference in the amount of heat required to increase the temperature of a sample and reference are measured as a function of temperature or time.

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- e-mail : cclock95@unist.ac.kr
- Phone : 052-217-4033

Basic Principle

Dynamic Mechanical Analysis measures the mechanical modulus of the material as function of temperature and time and frequency by giving the external force which vibrates to the specimen as stress of sinusoidal wave.

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- Autospec Premier (Waters, USA)

Specifications

- Sensitivity : 5×10^{-7} C/ μ g for methyl stearate at m/z 298.3, EI, 1,000 resolution and 5×10^{-8} C/ μ g for methyl stearate at m/z 298.3, EI, 10,000 resolution
- Mass range : 2-1,200 Daltons or more at full sensitivity
- Resolution : Continuously variable above 60,000 (10% valley) or more.
- Scan rate : Better than 3 scan/second over the calibrated mass range 500-50-500 or equivalent
- Ion source : EI plug in ion volume
- Mass Analyser Ion Optics : EBE Tri-sector geometry
- Vacuum system : Diffusion or Turbo pump.

Manager : Kim Chul-Su

- e-mail : kimcs@unist.ac.kr
- Phone : 052-217-4067

Basic Principle

Autospec premier is the same system being GC/MS, and a double-focusing magnetic sector mass spectrometer using Waters Micromass patented tri-sector (EBE) double focusing geometry and wide gap magnet providing high sensitivity, high resolution, and low background noise. The EBE design has an electric sector on each side of the magnetic sector. The first electric sector has demagnifying optics, which give high resolution and sensitivity. The second electric sector improves abundance sensitivity automatically rejects metastable ion interferences, and reduces background noise..

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Inductively Coupled Plasma / Mass Spectrometer

Model

- ELAN DRC-e (Perkin Elmer, USA)

Specifications

- Dual turbo-molecular pumps with ceramic bearings
- Vacuum levels: 1×10^{-6} Torr
- 27.12MHz ~ 40MHz RF Generator
- Dynamic Reaction Cell
- Thermal stabilization of the quadrupole power supply (QPS) at >2.5MHz Mass Spectrometer

Manager : Son Hee-Sik

- e-mail : redson7598@unist.ac.kr
- Phone : 052-217-4068

GC / MSD

Model

- 7890A GC & 5975C MSD (Agilent, USA)

Specifications

- Ionization mode
 - Electron ionization
 - Positive/negative chemical ionization
- Mass range : 1.6 ~ 1,050 amu
- Mass scan speed : Max. 12,500 amu/sec
- Detector: Triple-axis HED-Electron Multiplier
- Vacuum pumping system : 262L/sec
- Library : Wiley 8th with NIST 2008 MS library

Manager : Kim Chul-Su

- e-mail : kimcs@unist.ac.kr
- Phone : 052-217-4067

Basic Principle

ICP/MS is mass spectrometer that use inductively coupled plasma by Ion source, and it makes high ionization efficiency by use high temperature plasma and offers mass spectrums. Also It has low detection limit and both outstanding reproducibility and precision.

ICP/MS system is component of ion source, mass detector, and interface that is made efficiently to inject samples in the instrument, in addition to vacuum, controller, and injector.

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Basic Principle

GC/MSD makes an effective combination for chemical analysis. Mixtures of volatile substances are separated by gas chromatography. Via a transfer line the fractions of each peak are transferred into the ion source where the ionization takes place. After acceleration, the ions pass through a magnetic-field which causes deviation of their linear flight curve according to their respective mass. After this magnetic field the detector counts the ions.

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Machine Shop

Lab	Equipment	Model	Maker	Page
CNC Machine Tool	CNC5-axis Machining Center	C40U	Hermle, Germany	64
	CNC3-axis Machining Center	B300V	Hermle, Germany	65
	CNC Lathe	TSL-6	S&T, Korea	66
	CNC Surface Grinding	DGS-630CNC	Daesan, Korea	67
General purpose Machine Tool	Lathe	TIPL-410	S&T, Korea	68
	Milling	STM-2VM	Stonic, Korea	69
	Band Saw	KDBS-200	Kyoungdong, Korea	70
	Electro Discharge Machine	ZNC200M	Jooje, Korea	71
Measuring Instrument	Coordinate Measuring Machine	PGS	Dukin, Korea	72
	Three-dimensional Measurement	NV-3000	Nanosystem, Korea	73
	Semi auto formtracer system	525-421K-1	Mitutoyo, Japan	74
	Multi-Component Dynamometer	2825A	Kistler, UK	75
	Laser interferometer	XL-80	Renishaw, UK	76
Powerful Microscope	MF-1010B	Mitutoyo, Japan	77	

E-beam Lithography

Model

- NB3 (NBL, UK)


Specifications

- Theoretical beam size : 2.1 nm @100 KeV, 7 nA
- Line width : < 5 nm
- Deflection : vector scan, 55 MHz
- Address grid resolution : 1 nm, 1 mm main field
- Beam voltage : 30 keV ~ 100 keV
- Writing area : 195 mm × 195 mm
- Substrate size : 5 mm ~ 200 mm
- Automation : 10 chucks automatic loading
- Repeatable : 20 nm over wafer

Manager : Lim Dong-Ju

- e-mail : djlim@unist.ac.kr
- Phone : 052-217-4032

Mask Aligner

Model

- MA6/BA6 (Suss Microtec, Germany)

Specifications

- High Efficiency & High Accurate UV Exposure System
- UV Lamp: Hg 350W
- UV400: 350~450nm (I-, H-, G-line)
- Gap adjustment accuracy: 1 μm
- Methods : Top Side Alignment / Bottom Side Alignment
- Alignment accuracy: 0.5 μm with TSA / 1 μm with BSA

Manager : Shim Young-Dae

- e-mail : syd0620@unist.ac.kr
- Phone : 052-217-4064

Basic Principle

As the primary electron beam penetrates the resist, elastic collisions spread the beam slightly. Back-scattered electrons from the substrate may be spread over several microns. The secondary electrons generated in this process impart energy to resist through inelastic collisions. This causes the physical and chemical changes in the resist that renders in more soluble than the unexposed resist.

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Basic Principle

Performs alignment and exposure by UV Lamp on the wafers and substrates coated by photo resist through mask with alignment keys

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- e-mail : syd0620@unist.ac.kr
- Phone : 052-217-4064

Basic Principle

Nanoimprinting is a cutting-edge technology enabling sub 100nm scale features to be produced less costly in a parallel way. Its fundamental principle is to transfer the master patterns defined in the stamp (or termed as mask) to deformable materials such as photoresist spun on substrate by mechanical press. Positive reliefs on the stamp deform the resist upon imprinting to form physical matches.

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Basic Principle

The device which deposits the uniformly thick photoresists by using high speed.

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Applications

- Adhesion Treatment on Substrates
 - 표면 점착력 처리
- Photo Resist Coating
 - Photo Resist 도포
- Soft / Hard Baking
 - Photo Resist의 Hard/Soft Baking

Dielectric RIE

Model

- LABSTAR-R (TTL, Korea)

Specifications

- SiO₂ etch
 - Process gases : CHF₃, Ar
 - Clean gases : SF₆, Ar
 - Etch rate : > 35nm /min
 - Uniformity : (max-min) / (2X average) ±5% across 200mm diameter
 - Selectivity : SiO₂ : Resist = 3~5 : 1 (Resist dependent)
- SiN etch
 - Process gases : CHF₃, O₂
 - Clean gases : SF₆, Ar
 - Etch rate : > 50nm /min
 - Uniformity : (max-min)/(2X average) ±5% across 200mm diameter
 - Selectivity : SiN : Resist = 1.5~3 : 1 (Resist dependent)

Manager : Shim Young-Dae

- e-mail : syd0620@unist.ac.kr
- Phone : 052-217-4064

Metal RIE

Model

- LABSTAR-R (TTL, Korea)

Specifications

- Reactor(Process chamber) module :
 - 200mm electrode with water cooled by (heat exchange - temp. range: 5 ~150 C)
 - Capacity of sample holding: 4 and 6 inch wafer
- Primary Vacuum Loadlock module :
 - Including end-effector for transferring a wafer between loadlock and process chamber without leaking vacuum.
- MFCs with bypass lines (BCI3 / Cl2 / Ar / O₂ / SF₆) installed
- RF Generator & auto match network
 - 600W, 13.56MHz solid state RF Generator with automatic match unit for bias(RIE) power (Should use one of AE, ENI or Seren)
- Pt etch process : Etch Rate(10 nm /min or higher)
 - Selectivity to PR mask : Better than 1.5:1
 - WIW,WTW Uniformity : ±5%
- Au etch process : Etch Rate(25 nm/min or higher)
 - Selectivity to PR mask : Better than 2:1
 - WIW,WTW Uniformity : ±5%

Manager : Shim Young-Dae

- e-mail : syd0620@unist.ac.kr
- Phone : 052-217-4064

Basic Principle

Reactive Ion Etcher which forms fine patterns by dry etching placing wafer inside the process module and generating plasma.

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- Si / SiO₂ / Si₃N₄ Dry Etching
Si / SiO₂ / Si₃N₄ 건식 식각

Basic Principle

This system designed for various metal films, Al/Ti/Cr/W and Au etch process with corrosive chemistries.

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- V15-G (KAMI, Germany)

Specifications

- Dimension: W600 X H1,850 X D880 mm
- 250 x 250 x 250mm, aluminum chamber
- Reactor with 200mm water cooled chuck
- Microwave power generator : 2.45GHz, from 100 up to max. 600W
- Gas inlet : 2-gas inlet pipes (top side mounted)
- O₂, Ar gas control unit (MFC) : ≥ 200ml/min
- Process storage capability of 10 main process with 16 sub-process, rotation table 19"

Manager : Shim Young-Dae

- e-mail : syd0620@unist.ac.kr
- Phone : 052-217-4064

Wet Station

Model

- Wet Station (Donghun Tech. Korea)

Specifications

- Chemical supply : Manual
- Chemical Temp.: 50~120°C
- Chemical drain : BOTTLE Drain (AIR VALVE 구성)
- PR WET STATION(2SET)
 - Ultrasonic BATH, DEVELOP BATH, Q.D.R Bath
- ACID WET STATION
 - SC-1,2 BATH, SPM BATH, BOE BATH, DHF BATH, H₃PO₄ BATH, Q.D.R Bath
- SOLVENT WET STATION
- ETCHING WET STATION
- MEMS WET STATION
- CMOS WET STATION

Manager : Shim Young-Dae

- e-mail : syd0620@unist.ac.kr
- Phone : 052-217-4064

Basic Principle

The equipment uses gas plasma selectively Coated Photo Resist for the semiconductor process. And plasma etching at the semiconductor fabrication Process or Polymer which occurs during process.

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Basic Principle

This device allows 6" wafer to progress into photore sist develop, remove and wet etching process after cleaning using DI water.

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- KHD-306 (KSM, Korea)

Specifications

- Heater Spec(3 zone) : 600°C~1200°C
(PROCESS 온도 800°C~1150°C)
- Wet Oxidation Thickness uniformity
 - 2000Å target, edge 5mm, 9 point
 - within wafer : 3%, wafer to wafer : 3%, Run to Run : 3%
- N+ Annealing(Dry Oxidation) Thickness uniformity
 - 1000Å target, edge 5mm, 9 point
 - within wafer : 3%, wafer to wafer : 3%, Run to Run : 3%
- P+ Annealing(Dry Oxidation) Thickness uniformity
 - 1000Å target, edge 5mm, 9 point
 - within wafer : 3%, wafer to wafer : 3%, Run to Run : 3%
- POCl₃ doping(phosphorous doping) Resistivity uniformity
 - within wafer : 5%, wafer to wafer : 5%, Run to Run : 5%

Manager : Kim Hyung-IL

- e-mail : hikim78@unist.ac.kr
- Phone : 052-217-4065

PECVD

Model

- PEH-600 (SORONA, Korea)

Specifications

- High Plasma Density, High Dep. rate PE-CVD
- Wafer Size : 6", Max Substrate Temperature : 400 °C
- Power Supply : 600 W, 13.56 MHz RF Generator
- Process Gas : N₂O, N₂, CF₂/O₂, pure SiH₄/NH₃
- Process Guarantee : 1000Å ~ 10um
- With in Wafer Uniformity : ±3%
- Wafer to Wafer Uniformity : ±3 %
- Uniform process gas flow distribution
- Robotic arm for wafer handling
- RF generator connected to upper electrode through a close coupled automatic matching network

Manager : Kim Hyung-IL

- e-mail : hikim78@unist.ac.kr
- Phone : 052-217-4065

Basic Principle

- Wet Oxidation : Annealing process with vapor in high temperature Chamber
- Annealing(Dry Oxidation) : Annealing process in high temperature Chamber.
- POCl₃ doping(phosphorous doping) : Impurity injection process for PoCl₃ doping in high temperature chamber

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Basic Principle


PEH-600 PECVD System is a multipurpose tool capable of deposition silicon oxide, silicon nitride, amorphous silicon using electron energy (plasma) as the activation method to enable deposition at a low temperature and at a reasonable rate. Machine is controlled by a PC that runs the PEH-600 software. This equipment can hold a 6" or 4" wafer for processing.

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- SRN-120D (SORONA, Korea)

Specifications

- Substrate Size : 4~6 inch Silicon Wafer
- Substrate Rotation Speed: 0 ~ 60 rpm
- RF Pre-cleaning : 300W, automatic RF matching
- Substrate Heating : 600°C \pm 3% (Wafer temp)
- Target : ITO, ZnO, SiO₂, Al₂O₃
- Target size : 4inch
- Number of Target : Max. 4 different targets
- Available multi-layer processing
- Power source : 13.56Mhz RF power processing
- High vacuum pumping using turbo pump
- System control PC
- Film uniformity : less than $\pm 5\%$ using oxide material

Manager : Kim Hyung-IL

- e-mail : hikim78@unist.ac.kr
- Phone : 052-217-4065

Basic Principle

UHV CVD System is the device which creates new qualified wafer based on chemical reaction between Gas and Sample by spraying a Gas to keep Base pressure 1×10^{-10} . UHV CVD System can make sample which user is expected to control the flow of Gas, the temperature of sample and vacuum pressure. Wafer and Gas of UHV CVD System will be replaceable, and can be used in any aspect of work procedure. The equal thickness of the thin film on UHV CVD system is the key feature of strength compared with other Deposition has. LEED on the UHV CVD system is able to analyze the nature of sample such like an amplification, pattern and etc (so on?) of sample which created from the system.

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Basic Principle

SRN-120D Sputtering system is a multipurpose tool generating plasma inside chamber by supply 13.56Mhz rf power to cathode and capable of deposition thin film of insulator on wafer by sputtering a insulator target by self bias voltage. This device was designed to enable multi layer processing and co-processing with a maximum four targets.

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- SRN-120M ((SORONA, Korea)

Specifications

- Substrate Size: 4~6 inch Silicon Wafer
- Substrate Rotation Speed: 0 ~ 60 rpm
- RF Pre-cleaning : 300W, automatic RF matching
- Substrate Heating : 300°C ± 3% (Wafer temp)
- Target : Ti/Ta/Cr/Cu/Ni/Al/Mo Metal Depositon
- Target size : 4inch
- Number of Target : Max. 4 different targets
- Available multi-layer processing
- Power source : 3kw DC power processing
- High vacuum pumping using cryo pump
- System control PC
- Film uniformity : less than ±5% using oxide material

Manager : Kim Hyung-IL

- e-mail : hikim78@unist.ac.kr
- Phone : 052-217-4065

Basic Principle

SRN-120M Sputtering system is a multipurpose tool generating plasma inside chamber by supply 3kw DC power to cathode and capable of deposition thin film of metal materials on wafer by sputtering a metal target by DC negative voltage. This device was designed to enable multi layer processing and co-processing with a maximum four targets.

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- WC-4000 (Woosung Hi-vac, Korea)

Specifications

- E-Beam source assembly (Au, Cr, Ni, Ti, Al, Cu, Pt, Pd)
 - 6cc, 6 pocket, 10" cf mounting
 - Power supply 6kw, 208v
 - Source control module : Remote package type
 - Low profile sensor package 2.75"
 - Programble sweep controller
- Thickness control
 - Multi film thickness control 2Ch
 - Frequency resolution : 0.5Hz@6.0MHz
 - 10film storage
 - SO-100 oscillator

Manager : Kim Hyung-IL

- e-mail : hikim78@unist.ac.kr
- Phone : 052-217-4065

Basic Principle

This equipment is thermal evaporating or e-beam deposition method equipment. Electron beam source is shown by the power supply to the hot filament electron beam by the magnetic field induced by magnet, if deposition materials are deposited by the collision of an intensive electronic material is heated to evaporate. Growth of thin films appears on the substrate. Using this kind of surface treatment and coating and can have many applications, a broad range of equipment. Thermal evaporator to the organic and metal deposition at the same time is the equipment.

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SAM (Self assembly mono-layer) Coater

Model

- AVC-150M (SORONA, Korea)

Specifications

- Superior Monolayer Properties
 - Contact Angle : 5 ° - 110 °(water)
 - Work of Adhesion : 3-100 uJ/m²
 - RMS Roughness : <0.2 nm
 - Thickness : 1 monolayer
 - Uniformity : ±3%
 - Repeatability : ±3%(wfr-wfr), ±3%(run-run)
- System specifications
 - Wafer Sample size : ≤ 300mm(depends on application)
 - Process Mode : Auto, Semi-Auto, Manual
 - Base Pressure : < 30 mTorr
 - Precursor delivery : Liquid direct injection lines w/precision meter(Patent pending)
 - Precursor Usage : <0.1 cc/batch (depends on application and throughput)
 - Surface Treatment : RF Plasma (100~300W)
 - Reactor cleaning by O₂ plasma
 - Gas Lines : O₂ and purge N₂

Manager : Kim Hyung-IL

- e-mail : hikim78@unist.ac.kr
- Phone : 052-217-4065

Basic Principle

AVC-150M can grow polymer mono-layer with very solid bonding force by forming O-H layer by natural oxygen on sample and supplied polymer source from vaporizer. Therefore, this polymer mono-layer by Vapor SAM equipment, reduces surface energy of samples and changes the surface characteristic from hydrophilic to hydrophobic. Also, The very uniformly thin grown polymer mono-layer, roughly about 1nm thickness, can realize uniform growth even for pattern with circuit line width of tens of nm.

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Applications

- Nano patterning of bio devices by imprinting
임프린트 바이오소자 나노패터닝
- Nano patterning of display panel by imprinting
임프린트 FPD 나노패터닝
- Nano patterning of polymer memory device by imprinting
임프린트 메모리소자 나노패터닝
- Nano patterning of information memory devices by imprinting

Parylene Coater

Model

- Parylene Coater (Alpha plus, Korea)

Specifications

- Dimer type : C
- Deposition Rate : Over 1-2um/hr
- Substrate Size : Max. Dia. 200mm
- Process temperature : Furnace(R/T to 1000°C±3°C) / Vaporizer(R/T to 300°C ±3°C)
- Thickness variation : Within ±10% (Center to Edge)
- Substrate loading : Manual

Manager : Kim Hyung-IL

- e-mail : hikim78@unist.ac.kr
- Phone : 052-217-4065

Basic Principle

The Dimer is inserted into the vaporizer in the powder form and evaporated between 120-175 °C Furnace chamber heated to approximately 650-690 degrees through the monomer is converted to the substrate. After the process, the temperature is recovered to room temperature creating a Poly-para-Xylylene film coating.

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- e-mail : syd0620@unist.ac.kr
- Phone : 052-217-4064

Basic Principle

Versatile functionality and analysis is only of use if high-quality original images can be obtained. The VHX has succeeded in producing sharper, super high resolution observation images by applying sophisticated technologies to every process from lighting to image generation. The VHX-600 not only provides advanced functions that enable ultra-deep and high definition observation, it can be easily operated by anyone. The VHX also meets a variety of user's requests for evaluation time reduction and quality improvement. Observation, recording and measurement can all be performed within the VHX-600.

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- P-6 (KLA_Tencor, USA)

Specifications

- Features & Configuration
 - Single 2D & 3D Scanning Profiler
 - 150mm Diameter Vacuum & Motorized sample stage or bigger
 - Manual or Auto Theta Sample Stage
 - Zoom Optic : Top View 185~750 or 115~465X
 - Automatic Step Detection and Multiple Cursor Positioning
- Specification
 - Scan Length : 150 mm
 - CCTV Zoom : Top View 185~750 or 115~465x
 - Scan Method : Bi-directional stage moving for scan
 - Scan Speed : 2um ~ 25mm / sec
 - Sampling Rate : 50, 100, 200 Hz
 - Vertical Range / Resolution : 13 um / 0.1Angstrom.
 - Stylus Force :Adjustable between 1~15mg
 - L-Stylus : 2um Radius 60 Degree

Manager : Shim Young-Dae

- e-mail : syd0620@unist.ac.kr
- Phone : 052-217-4064

Basic Principle

The P-6 stylus profiler and surface analysis system offers a combination of advanced features for process development and manufacturing control of scientific research, photovoltaic solar manufacturing, data storage, MEMS, opto-electronics, and other industrial metrology applications. The P-6 stylus profiler and surface analysis system utilizes numerous platform and measurement technologies from KLA-Tencor's most advanced semiconductor profiler systems for surface topography analysis - programmable scan stage, low noise, and high quality, high resolution long scans but in a smaller, more economical design for substrates up to 150mm.

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- ST4000-DLX (KMAC, Korea)

Specifications

- Wavelength Range : 400nm~800nm
- Thickness Measurement Range: 100Å~35μm
- Lens Turret: M5x, M10x
- Measuring Spot Size: 40μm, 20μm
- Size of sample stage: 200mm x 200mm
- Manual Focusing Module: Coaxial coarse & fine focusing controls

Manager : Kim Hyung-IL

- e-mail : hikim78@unist.ac.kr
- Phone : 052-217-4065

Basic Principle

Some portion of the light that is incident on the thin film surface is reflected from the upper surface of the sample, other portion is reflected from the interfaces between the thin film and the substrate. The light may be reflected from the interfaces between the thin films in case of multiple thin-film layers. As these lights are coherent, generated from the same light source, they interfere with each other and show constructive and destructive interference phenomena according to wavelength. And, these measured reflected lights show specific spectrum shapes in the wavelength region such as, sine patterns or valley-like, depending on thickness and optical constants. From these spectrum shapes, optical characteristics such as thickness, a refractive index and an extinction coefficient – can be obtained by assuming some appropriate values and adapting them systematically.

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Applications

- Semiconductor
- Dielectric materials
- Polymer
- LCD
- Optical Coating
- Recordable Material
- Etc.

4-Point Probe System

Model

- CMT-SR2000N (AIT, Korea)

Specifications

- Measurement range : 1 mohm/sq ~ 2Mohm/sq
- Measurement accuracy : ±0.5 % (Standard resistor)
- Measure pattern : ASTM, SEMI, Customer design
- Data analysis : Data map, Contour & 3D mapping
- Measurement unit : ohm, ohm/sq, ohm/cm(두께입력)
- Current Source : 10nA ~ 100mA
- Voltage : 0 ~ 2,000mV
- Substrate Size : 200mm(wafer)/ 140X140mm(square)
- 4 Point probe head : JANDEL probe

Manager : Kim Hyung-IL

- e-mail : hikim78@unist.ac.kr
- Phone : 052-217-4065

Basic Principle

Full automatic system to measure sheet resistance and resistivity silicon wafer, metal, and solar cell, etc. This system can be operated by itself, furthermore, perfect remote control is available using a PC and exclusive software, and it gives various data analyses

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- AR06DM (Aaron, Korea)

Specifications

- Substrate size : 2" to 6"
- Cutting materials : Silicon, glass, sapphire
- X-axis(chuck table horizontal movement)
- Work-piece width setting range: 0.01~160mm
- Cutting range: 160 mm or more
- Cut speed : 0.05~300 mm/s or more

Manager : Kim Hyung-IL

- e-mail : hikim78@unist.ac.kr
- Phone : 052-217-4065

Basic Principle

Performs fine alignment and permanent bonding between wafer and wafer by heat and pressure in a vacuum chamber. It can be used for various types of semiconductor process, especially for MEMS sensors and LED devices.

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Basic Principle

A device to cut Silicon wafer into small chip by highspeed rotating diamond blade.

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- C40U (Hermle, Germany)

Specifications

- Table
 - Size : over than $\Phi 800\text{mm}$
 - swivel Range: $+25^\circ/-110^\circ$ degrees or more
 - Max. Load capacity: over than 1,200kg
 - Speed-swivelling axis A: over than 25 rpm
 - Speed-rotary axis C : over than 30 rpm
- Work area
 - Traverse X-Y-Z: over than 850X700X500mm
 - Linear rapid traverse X-Y-Z: over than 45m/min
- Main Spindle
 - Speed : over than 18,000rpm
- Automatic Tool Changer
 - Magazine pockets: over than 38 pockets
 - Maximum tool length: over than 300mm
 - Maximum tool diameter: over than $\Phi 90\text{mm}$

Manager : Kwon Kang-Uk

- e-mail : kku1050@unist.ac.kr
- Phone : 052-217-4066

Basic Principle

A Numerically controlled milling machine, that, when the automatic tool changer is attached, has improved accuracy and efficiency. It is capable of performing up to 5-axis simultaneous machining.

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- B300V (Hermle, Germany)

Specifications

- Table
 - Size : over than 1,000X560mm
 - Max. Load capacity: over than 700kg
- Work area
 - Traverse X-Y-Z: over than 700X550X500mm
 - Linear rapid traverse X-Y-Z: over than 30m/min
- Main Spindle
 - Speed : over than 15,000rpm
- Automatic Tool Changer
 - Magazine pockets: over than 24 pockets
 - Maximum tool length: over than 300mm
 - Maximum tool diameter: over than $\Phi 80\text{mm}$

Manager : Cha Jae-Hoon

- e-mail : cjh614@unist.ac.kr
- Phone : 052-217-4069

Basic Principle

The CNC General milling machine and CNC milling machine equipped with an automatic tool changer enables a series of multiple processing tasks automatically.

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- DGS-630CNC (Dae San, Korea)

Specifications

- Capacity
 - Working Surface of Table : 650X300mm
 - Maximum Travel: 750X340mm
 - Distance from Table Surface to Wheel Center : 550mm
 - Electro-Magnetic Chuck: 600X300X110mm
 - Graduation of Hand Wheel: 0.001, 0.01, 0.1mm
- Grinding Wheel
 - speed : over than 1800rpm
 - Diameter X Width X Bore: $\Phi 305 \times 38 \times 127$ mm
- Motor
 - Grinding Wheel Spindle: 3.75kw/4p
- Weight : 3,100kg

Manager : Kwon Kang-Uk

- e-mail : kku1050@unist.ac.kr
- Phone : 052-217-4066

Basic Principle

The iteration of the servo motor rotation ball screw line by work piece movement is to convert into (ATC) redirect. The relative movement is intended, ultimately, to cut up the shape of the parts that are to be processed.

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Basic Principle

This Tool consists of a cylindrical grinding wheel instead rotating at high speed to work pieces into flat or grooved materials with a high level of precision.

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- TIPL-410 (S&T, Korea)

Specifications

- Capacity
 - Swing over bed: $\Phi 410\text{mm}$
 - Swing over carriage: $\Phi 220\text{mm}$
 - Distance between centers : 1,060mm
- Spindle
 - Spindle nose: KS B4022-A-6
 - Spindle taper: MT NO.6
 - Spindle speed: 60 ~ 1,500rpm
- Feed
 - Longitudinal feed: 0.039 ~ 0.541mm/rev
 - Cross feed: 0.019 ~ 0.271mm/rev
 - Metric thread: 0.5 ~ 7mm/p

Manager : Cha Jae-Hoon

- e-mail : cjh614@unist.ac.kr
- Phone : 052-217-4069

Basic Principle

This machine is primarily designed to fashion cylinders from metal. The machine is also able to create conical, oval, curved and "screw" shapes and perform operations such as cutting and grinding.

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- e-mail : cjh614@unist.ac.kr
- Phone : 052-217-4069

Basic Principle

The vertical milling machine uses circumferential cutting edge milling cutters that are placed at regular intervals. Materials can be adjusted vertically, horizontally and back and forth. By moving the feed table.

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- KDBS-200 (KyongDong, Korea)

Specifications

- Capacity
 - Round bar: Φ 220mm
 - Rectangle: 220×400mm
 - Type: Vise Turning
 - Round bar: Φ 160mm
 - Recrangle: 165×225mm
- Saw Blade
 - Length: 3035mm
 - Width: 25mm
 - Thickness: 0.9mm
 - Blade Speed: 25, 40, 51, 71m/min
- Electric
 - Drive Motor : 1kW

Manager : Cha Jae-Hoon

- e-mail : cjh614@unist.ac.kr
- Phone : 052-217-4069

Basic Principle

This ring-shaped band saw rotates at high speed to cut metal, wood, etc.

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- e-mail : cjh614@unist.ac.kr
- Phone : 052-217-4069

Basic Principle


Electro discharge Machining (EDM) is based on the principle that. The amount of insulation from the continuous erosion of the spark discharge which is generated by the conductor (usually a metal work piece), can pierce a hole or a pre-programmed cut.

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- e-mail : kku1050@unist.ac.kr
- Phone : 052-217-4066

Basic Principle

This third dimension measuring instrument measures numerical values precisely, digitalizing an object's position in space. Through the use of a probe.

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- e-mail : kku1050@unist.ac.kr
- Phone : 052-217-4066

Basic Principle

The third dimension measuring instrument is measures numerical values precisely, digitalizing an object's position in space, with out probe coming into physical contactless with the object. Instead measure meals are assessed through the use of lenses and halogen light.

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- e-mail : kku1050@unist.ac.kr
- Phone : 052-217-4066

Basic Principle

The form tracer is a surface roughness tester, and surface outline Measuring device embedded in a single unit, thus saving cost and space, while providing effective tests. In total, two kinds of surface roughness and facets are measured.

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- 2825A (Kistler, UK)

Specifications

- Measuring range - Fx, Fy, Fz: -250~250N
-Mx, Mz : -11~11 N
- Overload (Fx, Fy, Fz): -300~300N
- Threshold: <0.002N
- Sensitivity - Fx, Fz: -26 pc/N - Fy: -13 pc/N
- Linearity, all ranges: <+/-0.4% FSO
- Hysteresis, all ranges: <+/-0.5 FSO
- Crosstalk: <+/- 2%
- Degree of protection EN60529: IP67
- Weight: - Dynamometer: 0.87Kg
-Top plate: 0.75Kg

Manager : Cha Jae-Hoon

- e-mail : cjh614@unist.ac.kr
- Phone : 052-217-4069

Basic Principle

This instrument measures the power of machinery. When the charge to the voltage amplifier is raised the signal transition occurs through the final measurement, which will be converted into a physical description.

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- XL-80 (Renishaw San, U.K)

Specifications

- Accuracy : ± 0.5 ppm
- Linear resolution : 1nm
- Maximum travel velocity : 4m/s
- Between each automatically updated environmental compensation: 7seconds
- Dynamic capture rate : 50kHz
- Linear range as standard: 80m

Manager : Kwon Kang-Uk

- e-mail : kku1050@unist.ac.kr
- Phone : 052-217-4066

Basic Principle

This machine measures accuracy of movement in other machines through the use of a laser (Helium and neon gas-combination). Due to the laser wavelength, precise measurement is possible.

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- MF-1010B (Mitutoyo, Japan)

Specifications

- Measuring range : X: 100mm, Y: 100mm
- Z-axis height of the maximum mass: 150mm
- Effective Class size: 180mm X 180mm
- Maximum weight of mass: 10Kg
- Digital counter resolution :
0.001/0.0005/0.0001mm
- Body(W X D X H) : 562 X 730 X 667mm

Manager : Kwon Kang-Uk

- e-mail : kku1050@unist.ac.kr
- Phone : 052-217-4066

Basic Principle

This microscope has a significant range of motion and allows objects to be viewed as in a standard microscope, also in micrometers.

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


How to get to UNIST

Location

- UNIST is located at Banyeon-ri in the western suburb of Ulsan. The city of Ulsan is located at 360 km away south of Seoul and 40 km away west of Pusan. The city can be easily reached by Seoul, Busan and other cities throughout the nation.

Address : 100 Banyeon-ri, Eonyang-eup, Ulju-gun, Ulsan Metropolitan City
울산광역시 울주군 언양읍 반연리 100

Airplane

- After arriving at Ulsan Airport, you can take a taxi to get to UNIST. The fare for taxi from Ulsan Airport to UNIST is around 15,000 won depending on the traffic situation.
Buses that run via the school are as follows : **No.133, 233, 433, 733, 807**

Train

- After arriving at Ulsan station, you can take a taxi to get to UNIST. The fare for taxi from Ulsan station to UNIST is around 15,000 won depending on the traffic situation.
Buses that run via the school are as follows : **No.133, 233, 433, 733, 807**

Express Bus

- After arriving at Ulsan Express Bus station, you can take a taxi to get to UNIST. The fare for taxi from Ulsan Bus Station to UNIST is around 15,000 won depending on the traffic situation.
Buses that run via the school are as follows : **No.133, 233, 433, 733, 807**

Car

- From Seoul, Daejeon, Daegu, Busan
Gyeongbu Expressway (Direction of Busan) → West Ulsan IC → Highway No.24 → UNIST

May I help you?


If you want to use our equipment,
we provide it reasonably.

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How to get to UNIST

Location

• UNIST is located at Banyeon-ri in the western suburb of Ulsan. The city of Ulsan is located at 360 km away south of Seoul and 40 km away west of Pusan. The city can be easily reached by Seoul, Busan and other cities throughout the nation.

Address : 100 Banyeon-ri, Eonyang-eup, Ulju-gun, Ulsan Metropolitan City
울산광역시 울주군 언양읍 반연리 100

Airplane

• After arriving at Ulsan Airport, you can take a taxi to get to UNIST. The fare for taxi from Ulsan Airport to UNIST is around 15,000 won depending on the traffic situation.
Buses that run via the school are as follows : **No.133, 233, 433, 733, 807**

Train

• After arriving at Ulsan station, you can take a taxi to get to UNIST. The fare for taxi from Ulsan station to UNIST is around 15,000 won depending on the traffic situation.
Buses that run via the school are as follows : **No.133, 233, 433, 733, 807**

Express Bus

• After arriving at Ulsan Express Bus station, you can take a taxi to get to UNIST. The fare for taxi from Ulsan Bus Station to UNIST is around 15,000 won depending on the traffic situation.
Buses that run via the school are as follows : **No.133, 233, 433, 733, 807**

Car

• From Seoul, Daejeon, Daegu, Busan
Gyeongbu Expressway (Direction of Busan) → West Ulsan IC → Highway No.24 → UNIST

UNIST
Ulsan National Institute of
Science and Technology

100 Banyeon-ri, Eonyang-eup, Ulju-gun, Ulsan Metropolitan City
Tel : 052-217-4032 ~ 4038, 4064 ~ 4069 / Fax : 052-217-4039
www.unist.ac.kr