HRXRD XRD

GID Manual

August, 2022

CONTACT

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

The screenshot shows the DIFFRAC plus XRD Commander software interface. The 'Requested' column for Theta, 2Theta, Phi, Chi, X, Y, and Z is highlighted with a red box. The 'Requested' values are 0 for Theta, 2Theta, Phi, Chi, X, and Y, and -1 for Z. The 'Move Drive' button is circled in red in the top menu bar.

Std	Aux	Actual	Requested
Theta	0	0	0
2Theta	0	0	0
Detector			
Div. Slit			
Antis. Slit	0.1	0.1	0.1
Phi	0	0	0
Chi	0	0	0
X	0	0	0
Y	0	0	0
Z	-1	-1	-1

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
The screenshot shows the DIFFRAC plus XRD Commander software interface. The 'Absorber' dropdown menu is highlighted with a red box and a circled '1', showing a list of options with '5283' selected. The 'Set' button next to it is also highlighted. The main plot area is empty, and the status bar at the bottom shows 'Drive Lock Disabled'.

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The screenshot displays the DIFFRAC plus XRD Commander software interface. The main window shows a diffraction pattern plot with a single sharp peak. The control panel at the bottom is highlighted with a red box and contains several numbered callouts:

- 1: Scantype dropdown menu (set to Detector Scan)
- 2: Start input field (set to -1)
- 3: Increment input field (set to 0.01)
- 4: Stop input field (set to 1)
- 5: Scanspeed input field (set to 0.1)
- 6: Sec / Step button
- 7: Start button (with waveform icon)
- 8: Cont button

The software interface also shows various parameters on the left side, including Theta, ZTheta, Detector, Div. Slit, Antis. Slit, Phi, Chi, X, Y, Z, Absorber, Rotation, Detector Slit, Shutter, and X-Ray settings.

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